ENG

FOREHEAD SUPPORT

The I.M.O. forehead support is an additional device that occasionally replaces the chin plate when this is removed to enable the patient to eat, shave, etc.

It consists of a polyethylene back plate with two press-studs that attach to the occipital support and an adjustable front band that when placed on the forehead prevents the head from moving forward. This band remains in position due to its Velcro closure.

PRELIMINARY OBSERVATIONS

- The doctor is the only person who can determine the duration of treatment.
- You must follow the specific instructions regarding the daily duration of use.
- The initial fitting should be conducted by your doctor or by an orthopaedic technician.
- The use of this product with topical medication must be supervised by your doctor.
- Where pressure is applied, the skin should not be damaged or hypersensitive.
- For hygiene and comfort it is advisable to wear a cotton garment underneath the product.
- Periodically check the condition of the belt and/or elastic straps to ensure they are tight and in good condition (free from grease, cuts, rips, etc.).

WASHING INSTRUCTIONS

Clean all surfaces with a cloth and soapy water to eliminate any dirt. Keep the Velcro® clean to maintain its adhesion characteristics by removing threads, lint, etc.

ASSEMBLY INSTRUCTIONS

- Fit chest plate A. Modify the shape if necessary by removing occipital plate B, disc C and chin support bar D. Bend horizontally to follow the contour of the chest and vertically fit it to the sternum.
- Reassemble disc C to chest plate, using the inner position for small patients and the outer position for larger patients.
- Adjust the angle of both shoulder supports so that each passes over the shoulder comfortably between the trapezius insertion and the acromion. Tighten the screws that hold the shoulder support discs until completely locked. Cross the straps between the scapulae and fasten them to the buckles of the chest plate.
- Fit chin plate D by positioning the pad under the chin in the correct position. Fix in position by placing the proper hole over the pin. Mark the selected hole so that the chin plate can be placed in the proper position again if it is removed. Close firmly.
- Fit occipital plate B by adjusting stops E at the proper height, inserting rod stops in the slots on the chest plate. Close the fasteners and bend posteriorly to position the occipital pad under the occiput.

CAUTION: Be careful when you shape the occipital bars. Do not deform in the area that joins to occipital plate B.

Note: Any serious incident involving the product must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.

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