EL

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ENG

Indications:

Post-operative shoulder treatments, post-rotator cuff and anteroposterior capsular reconstruction. Post shoulder surgery. Humeral head fractures. Post shoulder arthroscopy.

Precautions:

- Read these instructions carefully before using the orthosis with the accessory. If you have any questions, consult your doctor or the shop where you purchased it.
- The general instructions of use and specific indications determined by the doctor should always be followed.
- The orthosis requires prior assembly and fitting by an orthopaedic technician or qualified medical professional.
- For the product to perform its function, the right size for the patient should be used.
- For reasons of hygiene, this orthosis should only be used by one patient. The product should not be reused by other patients.
- The monitoring of treatment with this orthosis and its proper use should be reviewed frequently by your doctor, an orthopaedic technician or qualified medical professional.
- Your doctor should be informed if you experience any kind of skin condition, sensitisation or side effect.
- Before using topical creams with this orthosis, consult your doctor, an orthopaedic technician or qualified medical professional.
- In areas of support where pressure is applied, the skin should not be damaged or hypersensitive.
- For the orthosis to function properly, its components must be checked regularly by your doctor to ensure that they are in good condition.
- Check that there are no rips, tears or wear defects.
- Clean the product regularly and maintain proper personal hygiene.
- Make sure that the Velcro fasteners are free from fluff and lint to ensure that their closure properties are preserved.
- The medical staff responsible for monitoring treatment may explain the suitability of the product to the patient or suggest its replacement if it is damaged or worn.


Side effects:

Possible carotid hypersensitivity syndrome triggered by the use of slings in obese and/or short-necked patients. Avoid excessive pressure of the brace in the carotid sinus region.

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

Assembly instructions:


Components:

- 1: Bilateral abduction pillow 910 10°.
- 2: Waist belt.
- 3: Shoulder strap.
- 4: Sling.

Fitting the orthosis:

- Place the 910 10° abduction pillow (1) on the side of the body with the injured shoulder. Attach the waist belt (2) to the pillow (1), wrap around the body and fasten to the other side of the pillow with the Velcro® located on the end of the belt.
- Adjust the pillow (1) to the required height. Any excess waist belt (2) can be trimmed off.
- Attach the shoulder strap (3) next to the sling (4).
- For greater patient comfort, ensure that the pad is correctly positioned.

Product maintenance - washing instructions:

Make sure that the Velcro® fasteners are free from fluff and lint to ensure that their closure properties are preserved. Close all of the Velcro® fasteners before washing. Wash at or below 40°C. Do not bleach. Do not tumble dry. Do not iron.

Note: Any serious incident involving the product must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.

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