PSL150

Férula antebraquial

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## ENG

### **Description and intended use:**

Forearm splint composed of two sections joined together by adjustable straps with hook-and-loop fastening. Controls and protects the forearm preventing dorsal flexion movements of the wrist while allowing palm flexion. Interior padding with Plastazote® for enhanced comfort.

### **Indications:**

Treatment of forearm and wrist fractures. Postoperative.

### **Precautions:**

- Read these instructions carefully before using the splint. If you have any questions, consult your doctor or the shop where you purchased it.
- For the splint to be effective, it must be correctly fitted to the patient's morphology.
- Always follow these general instructions of use and any particular guidelines indicated by your general practitioner.
- Only a doctor is qualified to prescribe and determine the duration of treatment, and carry out subsequent monitoring.
- If you notice any side effects, skin conditions or sensitivity, remove the splint and consult your doctor immediately.
- Your doctor should consider the use of topical creams in conjunction with the splint.
- In areas where pressure is applied, the skin should not be damaged or hypersensitive.
- Although the orthosis is not for single use, it should only be used by one patient. The product should not be reused by other patients.
- The first fitting of the orthosis should be carried out either by qualified personnel (doctor, orthopaedic technician, etc.) or under their supervision.
- For the orthosis to function properly, its components must be regularly checked to ensure they are in good condition. The medical staff responsible for monitoring the treatment may recommend the product to the patient or suggest that it be replaced if it is damaged or worn.
- To ensure they function correctly, clean the hook-and-loop fastenings regularly, removing any material that may have stuck to them.
- Clean the product regularly and maintain proper hygiene.
- Do not wear the orthosis while having a bath, shower, etc.
- Do not let children play with this device.
- Some of the materials that make up the product are not flame retardant. Therefore, its use is not recommended in environments with excessive heat, fire or radiation.
- For minor irritations caused by sweat, we recommend using the cotton sleeve-type garment included in the package to separate the skin from contact with the splint.

As the product contains ferromagnetic components, the splint should not be used during magnetic resonance tests or radiation associated with diagnostic or therapeutic procedures.

### **Contraindications:**

Contraindicated for open scars with swelling, redness and accumulation of heat.

### **Instructions for proper fitting:**

To achieve maximum therapeutic effectiveness for different conditions and prolong the product's usable life, the appropriate size must be selected for each patient. As excessive compression can lead to intolerance, resulting in bruising or lack of colour in the limb, we recommend adjusting the level of compression so that it is firm but comfortable.

Any modifications to the product must be carried out by an orthopaedic technician or a legally trained health professional, and the end-user or person responsible for fitting the product must understand how it works and how to use it correctly.

- Fold the end of the fabric sleeve and cut an opening for the thumb, as shown in figures 1 and 2.
- Next, pull the fabric sleeve over the hand and forearm.
- Open the splint as shown in figure 4.
- Place the forearm on the splint, as shown in figure 5.
- Close the splint as shown in figure 6.
- Adjust the straps as shown in figure 7, and check that the orthosis is positioned as shown in figure 8.

**Note:** If the rigid sections need to be adjusted to fit the patient better, technicians can shape the outer part of each section (where there is no Plastazote®) using a heat gun at up to 180°C from a distance of 15 cm.

### **Recommendations-Warnings:**

This product is recommended for the treatment of specific conditions and is only suitable for the purposes indicated in these instructions or by your doctor.

### **Product maintenance - washing instructions:**

When the product is not in use, store it in its original packaging, in a dry place and at room temperature. Rigid part: Clean regularly using a damp cloth.

Fabric part: Hand wash with warm water (maximum 30°) and neutral soap. To dry the orthosis, use a towel to absorb as much moisture as possible and then leave to dry at room temperature. Do not iron or expose it to direct heat sources such as stoves, tumble dryers, direct sunlight, etc. When wearing or cleaning the product, do not use abrasive or corrosive substances, alcohols, ointments, or solvent fluids.

**Note:** Any serious incident involving the product must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.

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PRIM, S.A.