airtex^{PRO}

ATX730 · ATX721

Airtex rodillera rom	Airtex rodillera policéntrica
أير تكس مشدد علامة الر كبة روم	أير تكس مشدد علامة الر كبة متعدد المراكز
Airtex kolenní ortéza rom	Airtex kolenní ortéza polycentrická
Airtex kniebandage rom	Airtex polyzentrische kniebandage
Airtex knæbandage rom	Airtex polycentrisk knæbandage
Airtex επιγονατίδα rom	Airtex πολυκεντρική επιγονατίδα
Airtex rom põlveortoos	Airtex polütsentriline põlveortoos
Airtex rom knee support	Airtex polycentric knee support
Airtex-polvituki rom	Monenkeskinen airtex-polvituki
Airtex genouillère rom	Airtex genouillère polycentrique
Steznik za koljeno airtex rom	Policentrični steznik za koljeno airtex
Airtex rom térdrögzítő	Airtex policentrikus térdrögzítő
Airtex ginocchiera rom	Airtex ginocchiera policentrica
„Airtex“ kelio įtvarys rom	„Airtex“ policentrinis kelio įtvarys
AIRTEX celgala fiksators rom	AIRTEX policentrisks celgala fiksators
Airtex rom kniebrace	Airtex polycentrische kniebrace
Airtex rom-knebeskyttelse	Airtex polysentrisk knebeskyttelse
Airtex orteza stawu kolanowego typu rom	Airtex: orteza stawu kolanowego z przegubami policentrycznymi
Airtex joelheira rom	Airtex joelheira policêntrica
Airtex genunchieră rom	Airtex genunchieră policentrică
Airtex наколенник rom	Airtex наколенник полицентрический
Kolenná ortéza airtex s klbom rom	Polycentrická kolenná ortéza airtex
Airtex knäortos rom	Airtex policentrisk knäortos

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Description:

ATX730 - Semi-open tubular knee support with upper and lower tightening straps and patella pad, plus stabilisation and control of the knee's flexo-extension thanks to the joints with ROM system. Adjustable range of extension: 0°, 15°, 30°, 55°, 75°, 95° and 110°. Adjustable range of flex: 15° and 75°. Fits left and right.

The knee support comes with an Allen key and two pins.

ATX721 - Semi-open tubular knee support with upper and lower tightening straps, patella pad and bilateral aluminium polycentric joint to reinforce mediolateral stability.

The configuration of the tightening straps makes it possible to treat a number of different pathologies. Both straps can be crossed at the front or at the back, depending on the ailment being treated.

Indications - Intended use:

ATX730 - Control of knee mobility. Post-surgical treatment of the knee, collateral and patellofemoral ligament instability. Anteroposterior and mediolateral knee instability, anterior cruciate ligament (ACL) and posterior cruciate ligament (PCL) injury and post-injury rehabilitation.

ATX721 - Treatment of joint pain, bursitis, tendinitis, arthritis and osteoarthritis, patellar injuries, lateral knee ligament injuries, Osgood-Schlatter disease, chondromalacia patella, knee subluxation, post-injury and post-surgical treatment. Patellofemoral or lateral ligament instability, anterior cruciate ligament (ACL) and posterior cruciate ligament (PCL) injury.

For both models, all of the indications apply to any of the four fitting options, with the exception of treatment for ACL (option 1) and PCL (option 2).

Contraindications:

- Do not use on open wounds with swelling.
- Do not use if the area where the orthosis is placed is infected.
- Do not use in the case of skin allergies due to contact with biocompatible materials.
- Do not use if the area starts to redden and/or becomes excessively hot.

Warnings:

- Read these instructions carefully before using the orthosis. If you have any questions, consult your doctor or the shop that sold you the product.
- For the orthosis to serve its purpose, it must be correctly adjusted to the patient's morphology. The first fitting and adjustment must be carried out by qualified personnel (doctor, orthopaedic technician, etc.).
- Always follow these general instructions of use and any particular guidelines indicated by your general practitioner.
- Only a doctor is qualified to prescribe and determine the duration of treatment and carry out subsequent monitoring.

- If you notice any side effects, skin conditions or sensitivity, consult your doctor immediately.
- Your doctor might consider the use of topical creams in conjunction with the orthosis.
- In areas where pressure is applied, the skin should not be damaged or hypersensitive.
- Although the orthosis is not for single use, it should only be used by one patient. The product should not be reused by other patients.
- For the orthosis to function properly, its components must be regularly checked to ensure that they are in good condition. The medical staff responsible for monitoring the treatment may recommend the product to the patient or suggest that it be replaced if it is damaged or worn.
- To ensure they function correctly, clean the hook-and-loop fastenings regularly, removing any material that may have stuck to them.
- Clean the product regularly and maintain proper hygiene.
- Ask your doctor if you can remove the orthosis when bathing or showering. If you cannot, make sure you cover it or keep it away from the water.
- Do not let children play with this device.
- Do not change the position of the pins as established by the doctor.

Fitting instructions:

- It is essential that you select the correct size for each patient.
- You must make sure that the thigh area and all of the straps are open before starting to fit the orthosis.
- Place the orthosis on the leg so that the knee opening is over the kneecap, and so that the mechanical axis of the joints coincides with the anatomical axis of the knee and allows it to flex (Fig. 1)
- Then fasten the thigh section of the knee support (Fig. 2)
- When tightening, the straps must be crossed, whether at the front or at the back, in accordance with the pathology being treated. (Fig. 6-9)
- In order to be crossed at the back, the strap must be folded over itself (Fig. 5)
- The straps must be threaded through the guides on the opposite side to the one they are attached to. It may be necessary to reposition the direction of the strap's micro-hook fastener when crossing them. If necessary, the length of the strap can be adjusted (Fig. 10)

KNEE SUPPORT WITH FLEXO-EXTENSION ADJUSTMENT (ATX730):

Adjusting the ROM joint (ATX730):

You do not need to remove the joints from the knee support in order to adjust the angle limiter of the ROM joint. Using an Allen key, unscrew and remove the pin and insert it into the corresponding angle limiter. Once you have inserted the pin, screw it back in place.

If you want to limit flex and extension at the same time, insert the additional pin into the corresponding angle for each function.

Product maintenance - washing instructions:

This product is latex-free. Before washing, fasten all of the velcro fasteners. Wash with a neutral detergent at a maximum of 30°C. Do not dry clean. Do not iron. Do not use bleach. Do not tumble dry. Do not expose to direct heat sources such as stoves, heaters, radiators, direct sunlight, etc.

Remove the joints before washing, clean them with a damp cloth, then dry.

Note: Any serious incident involving the product must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.

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PRIM, S.A.