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ENG
Knee brace and intended use: Knee brace for the control of genu recurvatum consisting of a rigid 2.5 mm thick aluminium structure which is fixed to the femorotibial joint of the leg using length-adjustable anterior/ posterior straps which can be adjusted to the popliteal fossa. The structure is adjusted to the width of the leg with screws and allows a slight flexion of the knee during walking, while preventing hyperextension.

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AR

EN
Orthesis and intended use: Orthesis of knee for the control of genu recurvatum consisting of a rigid 2.5 mm thick aluminium structure which is fixed to the femorotibial joint of the leg using length-adjustable anterior/posterior straps which can be adjusted to the popliteal fossa. The structure is adjusted to the width of the leg with screws and allows a slight flexion of the knee during walking, while preventing hyperextension.

Indications: Brace for treatment of genu recurvatum. Stabilisation of the femorotibial joint in the anteroposterior plane, limiting hyperextension.

Warnings:

- Read these instructions carefully before using the orthosis. If you have any questions, consult your doctor or the shop where you purchased it.
- For the brace to perform its function, the right size for the patient must be used. The first fitting and adjustment of the brace must be carried out by qualified personnel (doctor, orthopaedic technician, etc.)
- Always follow these general instructions of use and any particular guidelines indicated by your general practitioner.
- Only a doctor is qualified to prescribe and determine the duration of treatment, and carry out subsequent monitoring.
- Your doctor might consider the use of topical creams in conjunction with the orthosis.
- Do not wear the orthosis while having a bath, shower, etc.
- Do not let children play with this device.
- Some materials of the product are not flame retardant, therefore it is recommended not to use it in environments with excessive heat, fire or radiation.

Fitting and adjustment instructions:

- The first step is to adjust the width of the brace and the tension of the popliteal girth.
1.1 Release the upper and lower front straps, unthreading one side for comfort. Loosen the screws at the back of the brace, without unscrewing them completely.
1.2 Place the brace over the patient's leg from behind, adjusting the central pad with the popliteal fossa Adjust the width and tighten the screws.
Important: The red circle in the joint area should be on the left and the blue circle on the right side of the knee.
1.3 Adjust the tension of the popliteal fossa strap by removing the pad and tightening the straps. Replace the pad.
2. Once the brace is adjusted, insert the leg by aligning the central pad with the popliteal fossa.
3. Adjust the front and rear straps to the desired tension by pulling on both sides of the straps so that all pads are centred. First adjust the rear straps, which remain fixed, and then adjust the front straps in the same way.
4. You can adjust the length of the front straps. Remove the micro snap hooks, use scissors to cut to the desired length and replace them.

Product maintenance – washing instructions:
To wash the textile components, remove the straps and pads from the metal frame and wash by hand with a neutral soap. Do not use bleach. Before washing, snap shut all micro-hooks; do not expose to direct heat sources such as stoves, heaters, radiators, direct sunlight, etc. Clean the metal structure with a damp cloth and then wipe dry.

Note: Any serious incident involving the product must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN

EN