mottio^{sport}

Cincha patelar

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1. Description and intended use:

The main function of the infrapatellar support is to relieve anterior knee pain using a pad that acts as a shock absorber, dispersing pressure and absorbing impacts on the patellar tendon. It has been designed to stabilise the kneecap, allowing the wearer to return to high-impact sports with greater protection and less discomfort. One size. It is recommended exclusively for patients with a knee circumference of between 25 and 46 cm.

2. Indications:

Treatment of Osgood-Schlatter disease. Treatment of patellar tendinitis (jumper's knee). Treatment of patellar chondromalacia.

3. Contraindications:

- Do not use on open wounds.
- Do not use if the area where the orthosis is placed is infected.
- Do not use in the case of skin allergies due to contact with biocompatible materials.
- Do not use if the area starts to redden and/or becomes excessively warm

- Do not use continuously throughout the day without breaks, as the lack of natural movement may lead to muscle weakness.
- Please consult your doctor before using these products if you suffer from poor circulation or nerve disorders.
- If you have any doubts about the suitability of the product, consult a doctor.

4. Warnings:

- Read these instructions carefully before using the orthosis. Always follow the general and particular instructions of use of your doctor. If you have any questions, consult your doctor or the supplier you purchased it from before use.
- Only a doctor is qualified to prescribe the orthosis and determine the duration of treatment. Treatment should be regularly monitored by qualified personnel.
- The orthosis requires assembly and prior adjustment by qualified personnel.
- Although the orthosis is not for single use, it should only be used by one patient. The product should not be reused by other patients. Clean the product regularly and maintain adequate personal hygiene.
- If you notice any side effects, skin conditions or sensitivity, consult your doctor immediately. Discontinue use until assessed by a specialist.
- Your doctor should consider the use of topical creams in conjunction with the orthosis.
- In areas of support where pressure is applied, the skin should not be damaged or hypersensitive.
- For the orthosis to function correctly, its components must be regularly checked to ensure they are in good condition. The medical personnel responsible for monitoring the treatment may recommend the product to the patient or suggest that it be replaced if it is damaged or worn.
- Ask your doctor if you can remove the orthosis when bathing or showering. If you cannot, make sure you cover it or keep it away from the water.

5. Instructions for use:

- Detach one of the parts of the adjustment strap, leaving it completely free. Fig. 1.
- Fit the strap to the leg at the height of the patellar tendon (just below the patella) and centre it.
- Pass the strap through the buckle again, pull both straps backwards around the leg and fasten. Fig. 2.
- Stand up and walk around. Check that the orthosis fits snugly and

comfortably, without causing any pain or discomfort. If it feels too loose or too tight, adjust the straps slightly.

6. Product care:

Before washing, close all the velcro fasteners. Hand wash with warm water and mild soap. Do not use bleach. Do not tumble dry. Do not iron. Do not dry clean. Do not expose the orthosis to direct heat sources such as stoves, heaters, radiators or direct sunlight.

Note: Any serious incident involving the product must be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located.

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