



Maintain your body heat. Breathe





neoprair



Maintain your body heat. Breathe



PRIM, S.A.
P. Industrial N.1
Calle F, 15
28938 Móstoles
(MADRID) • SPAIN

export@prim.es
www.prim.es



A dual effect that enhances well-being

Our new NEOPRAIR orthoses feature AIRLOCK SYSTEM technology to enable the skin to breathe while maintaining body heat. A dual effect that enhances comfort and pain relief.

The fabric, which provides multidirectional elasticity and three-dimensional compression, makes NEOPRAIR orthoses **100% lightweight** and adaptable to the shape and mobility of the joints, enhancing its ability to offer support and stabilisation during sporting activities.

NEOPRAIR. Heat and breathability.

SPORT SUPPORTS

100% lightweight



Maintains your body heat

Expels sweat

AirLock SYSTEM



Therapeutic heat

The internal heat cycle is maintained thanks to the intrinsic properties of neoprene.

Heat that relieves pain and reduces recovery time.

Breathability

Its AirLock System technology generates channels in the neoprene which enable sweat molecules to evaporate allowing the skin to breathe.

Maintain your body heat. Breathe





Closed knee support

- Inner silicone bands with anti-slip effect.
- Combines breathable neoprene of different thicknesses to provide greater comfort in the popliteal area.

Knee support with patella pad

- Patella pad to stabilise the kneecap.
- Inner silicone bands with anti-slip effect.
- Combines breathable neoprene of different thicknesses to provide greater comfort in the popliteal area.

Knee support with patella pad and splints

- Patella pad to stabilise the kneecap.
- 4 splints to improve lateral stability.
- Inner silicone bands with anti-slip effect.
- Combines breathable neoprene of different thicknesses to provide greater comfort in the popliteal area.

Knee support with splints and straps

- Patella pad to stabilise the kneecap.
- 4 splints to improve lateral stability.
- Double fastening strap for an optimum fit. Enables control of pressure.
- Inner silicone bands with anti-slip effect.
- Combines breathable neoprene of different thicknesses to provide greater comfort in the popliteal area.

Knee support with polycentric joint

- Patella pad to stabilise the kneecap.
- 2 polycentric aluminium joints to reinforce mediolateral stability.
- Double fastening strap for an optimum fit. Enables control of pressure.
- Inner silicone bands with anti-slip effect.
- Combines breathable neoprene of different thicknesses to provide greater comfort in the popliteal area.

Wraparound knee support with polycentric joint

- Open kneecap including a patella pad.
- 2 polycentric aluminium joints to reinforce mediolateral stability.
- Inner silicone bands with anti-slip effect.
- Opening in the popliteal area or greater comfort.
- Wraparound fastening system.



Lumbosacral back support with tensioner

- Secure front fastening.
- Removable tensioners for greater compression.
- 4 metallic stays to provide excellent containment. Removable.

NP139

S 35-38 / M 38-41
L 41-45 / XL 45-50
Bilateral



NP133

S 35-38 / M 38-41
L 41-45 / XL 45-50
Bilateral



NP137

S 35-38 / M 38-41
L 41-45 / XL 45-50
Bilateral



NP136

S 35-38 / M 38-41
L 41-45 / XL 45-50
Bilateral



NP140

S 35-38 / M 38-41
L 41-45 / XL 45-50
Bilateral



NP141

S 32-36 / M 36-41
L 41-49
Bilateral



NP156

S 80-90 / M 90-100
L 100-112 / XL 112-125



Cross-strap ankle support

- 3 tibial fastenings for an optimum fit.
- Double elastic band for figure-of-eight strapping.
- Combines breathable neoprene of different thicknesses to provide greater comfort in the plantar area. Open heel.

NP115

S 20-22 / M 22-24 / L 24-26
XL 26-29 / XXL 29-32
Bilateral



Figure-of-eight ankle support with stabilisers

- 3 tibial fastenings for an optimum fit.
- Double elastic band for figure-of-eight strapping.
- Lateral anatomical pieces for greater stabilisation without discomfort. Removable.
- Combines breathable neoprene of different thicknesses to provide greater comfort in the plantar area. Open heel.

NP114

S 20-22 / M 22-24 / L 24-26
XL 26-29 / XXL 29-32
Bilateral



Shoulder support

- Elastic strap with secure fastening for abduction adjustment.
- Armpit pad for greater comfort.
- Combines breathable neoprene of different thicknesses to provide greater comfort in the arm area.

NP174

S 26-29 / M 29-32
L 32-36
Left and right versions



Elbow support

- Forearm strap for selective adjustment.
- Combines breathable neoprene of different thicknesses to provide greater comfort in the elbow area.

NP173

S 23-26 / M 26-29
L 29-32 / XL 32-35
Bilateral



Wrist support with splints

- Removable aluminium palmar splint and plastic stays to enhance stabilisation.
- 4 secure fastenings for an optimum fit.
- Fastening strap on the wrist for pressure control.
- Free thumb.

NP162

S 15-17 / M 17-19
L 19-21 / XL 21-24
Left and right versions



Wrist and thumb support

- Removable aluminium thumb splint to help unload the metacarpal joints.
- 4 secure fastenings for an optimum fit.

NP163

S 15-17 / M 17-19
L 19-21 / XL 21-24
Bilateral



Heat and breathability