



*New sizes.
Better adaptation*

Humerux

Humerus orthosis




Maximum control

New sizes. Better adaptation

Humerux

Humerus orthosis



HUMERUX is a new orthosis that offers maximum control in humerus injuries, thanks to its innovative traction system.

Indicated for the treatment of diaphyseal humerus fractures, postoperative treatment of fractures and control of the instability of the fractures focus during the consolidation process.




Maximum control

Y-traction System uses a traction device to increase and adjust the tension, helping to maximise control over the alignment, fixation and immobilisation of the humerus.

Padded and breathable technical fabric: comfortable and respectful with the skin.

100% LATEX FREE



Valves and straps stabilize the humerus by applying uniform and personalized compression to the surrounding soft tissues.

Bivalve structure with removable, high density PE thermoformed plates: can be heat formed to achieve a perfect anatomical adaptation to the patient.

*The valves can be shaped with a heat gun at 150°C at about 15 cm. distance



Easy-fit chest strap to prevent the orthosis from sliding down. Padded armpit support for comfortable underarm fitting avoiding any type of friction or discomfort.

Fixation system using 3 humeral straps, threaded through 3 butterfly buckles with micro-hook closure and an anchor stop: they ensure an easy, comfortable and secure adjustment

Humerux EVOH300

Sizes

Short

S	18-26 cm	25 cm
M	26-34 cm	27 cm
L	34-42 cm	29 cm
XL	42-50 cm	32 cm

Sizes

Long

S	18-26 cm	30 cm
M	26-34 cm	32 cm
L	34-42 cm	34 cm
XL	42-50 cm	36 cm

A Range (upper arm contour)

B Length of the orthosis
Bilateral