




immobilising wrist supports

# airMed<sup>®</sup>

Cutting-edge functional immobilisation

Technologically advanced



High-frequency technology





immobilising wrist supports

# airMed®

**Cutting-edge functional immobilisation**

**AIRMED** is a range of wrist supports made using high-frequency technology and designed with clean, pure lines that contribute to the effectiveness of the treatment while maximising comfort.



**100% LATEX FREE** ✓

*Technologically advanced*

**Breathable, lightweight and adaptable fabrics**

*With perforated ventilation areas*

**Functional immobilisation splints**

*With anatomically contoured aluminium splints*

**"itself" fit system**

*Adjust the orthosis easily and safely using just one hand.*

High-frequency technology






## Immobilising thumb abduction orthosis

### AM202

Bilateral

	A Wrist circumference cm.	B Thumb circumference cm.
● AM202G S	14-18	5-7
● AM202G L	18-22	6-8
● AM202 S	14-18	5-7
● AM202 L	18-22	6-8

## Short wrist support

### AM203G Thumbless

	A Wrist circumference cm.	B Length cm.
● AM203G SR	14-18	16
● AM203G SL	14-18	16
● AM203G LR	18-22	16
● AM203G LL	18-22	16

Right and left versions (R,L)



### AM206G With thumb

	A Wrist circumference cm.	B Length cm.
● AM206G SR	14-18	16
● AM206G SL	14-18	16
● AM206G LR	18-22	16
● AM206G LL	18-22	16

Right and left versions (R,L)



## Long wrist support

### AM204G Thumbless

	A Wrist circumference cm.	B Length cm.
● AM204G SR	14-18	20
● AM204G SL	14-18	20
● AM204G LR	18-22	20
● AM204G LL	18-22	20

Right and left versions (R,L)



### AM205G Thumbless

	A Wrist circumference cm.	B Length cm.
● AM205G S	14-18	20
● AM205G L	18-22	20

Bilateral



### AM208G With thumb

	A Wrist circumference cm.	B Longitud cm.
● AM208G SR	14-18	20
● AM208G SL	14-18	20
● AM208G LR	18-22	20
● AM208G LL	18-22	20

Guarda lado (R,L)



## Thumb immobiliser accessory

### AM201G ● One Size

Can be used with  
AM203G, AM204G  
and AM205G

