

Multi-positional foot-ankle splint

PSL800



Controlled immobilisation

Anti-rotation splint
for foot-ankle positioning.
Immobilises the foot
and controls the posture
of the lower limb joints.

Indicated for postoperative therapy of the knee, hip and foot whenever control of lower limb rotation and dorsal flexion (neutral position, internal rotation or external rotation) is required. Hemiplegia, paraplegia, lower limb paralysis, immobilised or comatose patient (ICU patient). Plantar fasciitis.

Multi-positional foot-ankle splint PSL800

Controlled immobilisation



Adjustment system with double straps and buckles for closure. Integrated protective pads for optimum comfort.

Side splints. Allow plantar and dorsal flexion to be adjusted by up to 60°, in 10° intervals.



Comfortable padding. Three-layer, breathable and washable padded inner lining. Easy to remove and fit.



Plantar plate to stabilise the splint, with internal or external regulation from 0° to 15°.



100% LATEX FREE ✓



export@prim.es

Size and measurements

	Foot size	A	B
S	<32	28,5	26
M	32-38	30	28,5
L	39-44	32	30,5
XL	>44	34	32

Bilateral

Height (cm) A



B
Length (cm)