



Hip orthosis

PSL699 OS

Provides precise protection to the hip joint, stabilising and limiting its movement after surgery.

**Maximum stability,
Maximum control**



**articulation
full-control**

Hip orthosis

PSL699 OS



Maximum stability, Maximum control

Girdle-shaped pelvic band, with stabilising lumbar support and tensioners to adjust compression. Divided into 3 sections, in order to facilitate fitting for bed-bound patients.

Articulation *full Control*

- Enables flexion and extension movements to be adjusted in 15° increments, from 0° up to 90°.
- Controls abduction in the 0°, 15° and 30° positions.
- Offers 5 positions for personalising the angle of inclination of the upper bar.

Bivalve thigh cuff with condylar support

controls the leg's internal rotation and retraction by bringing it to external rotation.

2 telescopic plates

allow the orthosis to be adjusted to the height of each patient.

Removable protective padding

adds comfort and distributes pressure.





Hip orthosis

PSL699 OS

Maximum stability, maximum control

Indications: post hip surgery.
Replacement or repair of total hip prosthesis. Hip dislocation.
Fractures associated with hip replacement. Instability of the hip joint. Osteoporosis. Loosening of hip endoprotheses. Moderate to severe coxarthrosis.



breathable
latex-free

One size · bilateral

Waist circumference
range:

60 a 140 cm

Thigh circumference
range:

40 a 65 cm

