



Antebrachial orthosis

Forearm splint for functional treatment

PSL150

Bivalve orthosis

made in moulded thermoplastic, attached by straps with guide pins for a precise fit. *The upper part* is specially shaped to limit dorsiflexion of the wrist.



Inner lining

with Plastazote padding for added comfort.




100% cotton tubular

to prevent the injuries that can occur between the valves.



Adjustable straps

made in velour for an agreeable skin contact. Finished with micro-velcro and Velcro fasteners to enable safe and personalised adjustment.



1° Cut the strap



2° Incorporate to closure

Sizes and measurements



Size	Circumference	Length
S	20 - 24 cm	27,5 cm
M	24 - 28 cm	29 cm
L	28 - 32 cm	30,5 cm
XL	32 - 36 cm	32 cm

Right and left versions

100% LATEX FREE ✓