



PH 

A 

R 

M 

A

general catalogue





**PRIM**

**PHARMA**  
general catalogue





# Health, independence and well-being

Three vital aspirations always present in society.

And in pharmacies are found many of the keys to making them a reality.

**PRIM** would like to accompany you in the difficult and, at the same time, exciting work that you perform as a health professional.

# We have much



**+ 150**  
years of history

**PRIM Group** is a benchmark company in the Spanish market



Pharmacy



Orthopaedics



Physiotherapy and rehabilitation



Medical and hospital supplies



Spas



Listed on  
the stock  
market

Since 1985



**673**  
professionals

+14,8% in the last  
two years



**169,0M€**  
turnover 2021

+14,7% in the last  
five years



**2**  
own factories

We are a benchmark in  
orthotics and podiatry with o  
ur own manufactured products

# to offer you



**innovation**

**We anticipate the needs of the market** with the development of new products, innovating models and materials



**+1000 items**

**Your single supplier of** orthotics, podiatry, healthcare, technical aids and parapharmacy products



**281.247 orders**

**Managed in 2021**



**Global company**

**Export of own products** and search for new products for the local market across 5 continents




**Leaders in quality**

**Our products meet accepted manufacturing standards,** ensuring consistent quality



**8 regional offices**

**We have a well-established commercial structure with our own network** to provide a professional and individualised response to market needs



**Personal attention**

**We listen to your needs** to guarantee your business's present and open new opportunities for the future





# PRIM display furniture and accessories



Impact, visibility and sales.

# PRIM FLOOR STAND


Ref.: MBSUELOPRIMEXP  
C.I.: 828535

- Individual display stands to showcase product lines.

Can be joined together to create custom PRIM display areas: 2-3-4 stands.

Perforated structure allows hooks or trays to be placed at the desired height.

- 18 hooks included and additional bags of 4 hooks available on request (Ref.: MBGANCHOEXP).
- Trays not included: request separately (Ref.: MBBANDEJAEXP).
- Different signs available depending on the product line to be displayed: request separately.
- Connecting brackets not included: request separately depending on the number of displays to be connected (Ref.: MBUNIONEXP).



Floor stand  
AQTIVO SPORT



Floor stand  
AQTIVO SKIN



Floor stand  
AQTIVO MIX V



Floor stand  
COMFORSIL /  
COMFORGEL



Floor stand  
MOTTIO / AIRMED /  
AIRTEX /  
MYPRIM KIDS



Floor stand  
NEOPRAIR / MOTTIO /  
AIRMED / MYPRIM  
KIDS / AIRTEX



Floor stand  
PRIMSPINE



**CONNECTED FLOOR STAND**  
Ref.: MBUNIONEXP  
C.I.: 828536



**MANNEQUIN**  
Ref.: OTSCMANEXP  
C.I.: 820362



- Our iconic blue mannequin wears our product lines so the consumer can see them out of the packaging.

Complements the products displayed on the shelf.

Measurements: 200 x 50cm.  
(height x depth)

## PRIM WALL PANEL

Ref.: MBPANELPRIMEXP  
C.I.: 828537

- Individual panel display to hang on a wall.

Hang as many panels as the space at the point of sale allows: 2-3-4 panels.

Perforated structure allows hooks or trays to be placed at the desired height.

Brackets and wall template included.

- 18 hooks included and additional bags of 4 hooks available on request (Ref.: MBGANCHOEXP).
- Trays not included: request separately (Ref.: MBBANDEJAEXP).
- Different signs available depending on the product line to be displayed: request separately.



Wall panel display  
AQTIVO MIX H



Wall panel display  
NEOPRAIR / MOTTIO  
/ AIRMED / MY KIDS /  
AIRTEX



Wall panel display  
COMFORSIL /  
COMFORGEL



Wall panel display  
PRIMSPINE

## PRIM ISLAND DISPLAY

Ref.: MBISLAPRIMEXP

C.I.: 828538

- Free-standing island display stands to showcase products off-the-shelf in a prime position at the point of sale.

Cannot be joined together.  
Perforated structure allows hooks or trays to be placed at the desired height.  
36 hooks included.

- Trays not included: request separately (Ref.: MBBANDEJAEXP).
- Different signs available depending on the product line to be displayed: request separately.



Island display stand  
AQTIVO / PRIMSPINE

## PRIM DISPLAY STAND TRAY

Ref.: MBBANDEJAEXP  
C.I.: 828540



- Trays compatible with the different types of display stand, with attachment hooks.

Request separately depending on the requirements of the planned display.

Trays shipped separately



## LARGE FRAME

Ref.: MBMARCOLEXP  
C.I.: 828558



- They highlight and enhance framed products on the display, panel or island.

The 2 special hooks for placement (large or small frame) are supplied separately:

LONG FRAME HOOK  
Ref.: MBGANCHO2EXP  
C.I.: 828563



## SMALL FRAME

Ref.: MBMARCOSXP  
C.I.: 828559



## DISPLAY STAND SIGN



Display stand sign - PORT EXPORT  
Ref.: MBCARTELASPORTEX  
C.I.: 828543




Display stand sign - ORTHOPEDIC EXPORT  
Ref.: MBCARTELAORTOEX  
C.I.: 828544



Display stand sign - FOOT EXPORT  
Ref.: MBCARTELAFOOTEX  
C.I.: 828545

### • DISPLAY STAND SIGNS

Signs compatible with the different types of display stand: request separately for each product line to be displayed. (for island stands request 2 signs, one for each side of the display)





# 1

- 1.1. Myprim Kids
- 1.2. Orthotic treatment
- 1.3. Aqivo Sport
- 1.4. Aqivo Skin
- 1.5. Neoprair
- 1.6. Mottio
- 1.7. Airtex
- 1.8. Airmed
- 1.9. Essencial
- 1.10. PrimSpine Back Supports
- 1.11. Other orthoses

### New paediatric joint injury recovery line for children

A line with bright colours and fun shapes specially created for children, designed with the latest materials to ensure the comfort of the entire line, vital to the well-being of the child.



2/12  
years

breathable  
latex free



### Innovative development

MYPRIM KIDS has been developed with the latest bonding technology, which allows the seamless joining and adhesion of the different parts of the product.



### Absolute comfort

The absence of thread and seams significantly reduces chafing of the skin, preventing discomfort and increasing the well-being of the child.

## Display stands



DISPLAY STAND SIGN

PRIM FLOOR STAND -  
MBSUELOPRIMEXP

Floor stand  
MOTTIO / AIRMED /  
AIRTEX / MYPRIM KIDS



Shelf display -  
MPKREGLIN

Measurements  
40x9 cm.





move with freedom


## Children's supports



### Wraparound wrist brace - MPK100

Its flexible and padded design naturally wraps around the thumb and wrist. It exerts the necessary compression on the trapeziometacarpal joint to relieve pain while maintaining the functionality of the hand. It features velcro fasteners to regulate compression. Fits left and right.


Ref.	Size	Wrist circumference (cm)
MPK100 1	T1	10-12
MPK100 2	T2	12-15




### Immobilising wrist brace - MPK101

Immobilises the wrist joint thanks to an anatomically contoured palmar splint made from formed malleable aluminium, which limits palmar flexion. A dorsal splint provides the necessary stiffness to improve the stabilisation of the joint. Its 3-phase velcro fasteners enable the wrist brace to be easily and safely tightened. Left and right versions.

Ref.	Size	Wrist circumference (cm)
MPK101 1R	T1R	10-12
MPK101 1L	T1L	10-12
MPK101 2R	T2R	12-15
MPK101 2L	T2L	12-15




### Immobilising thumb attachment - MPK102

Immobilising thumb attachment to fit to the MPK101 wrist brace. Features a mouldable aluminium splint that enables the thumb joint to be immobilised as the injury requires. Fits left and right.


Ref.	Size
MPK102 OS	One Size



### Lumbosacral back support with tensioner - MPK600

Equipped with 2 steel stays at the back, it provides the necessary support for stabilising the spine and facilitating the unloading of the spinal structures. Features removable tensioners to provide greater stabilisation if required. Its velcro fastening enables the degree of compression to be regulated.


Ref.	Size	Waist circumference (cm)
MPK600 1	T1	50-60
MPK600 2	T2	60-70




### Clavicle immobiliser - MPK300

Immobilises and keeps the shoulders back thanks to adjustable non-elastic straps with buckle system. Covered with foam protection, they fit the armpit comfortably, preventing rubbing. Features an anatomically contoured cushioned dorsal pad for support and comfort. Its non-visible velcro fasteners can be removed to allow the length of the strap to be shortened.


Ref.	Size	Circumference under armpit (cm)
MPK300 1	T1	48-58
MPK300 2	T2	58-72




### Sling - MPK400

Made from a honeycomb technical sports material with soft terry cloth inside to absorb moisture and help maintain healthy skin. Features a cushioned humeral pad that minimises pressure on the neck and shoulder, providing great comfort. Adjustable strap and thumb loop. Ergonomic design. Fits left and right.

Ref.	Size	Forearm length (cm)
MPK400 1	T1	<21
MPK400 2	T2	21-25




### Shoulder immobilisation strap - MPK401

Shoulder immobilising attachment that fits to sling MPK400, achieving a 2 in 1. In the shape of a strap, it wraps around the arm and trunk, preventing the movement of the shoulder joint, if so required. Fits left and right.

Ref.	Size
MPK401 OS	One Size

## Children's supports



### Pavlik harness Child hip abduction splint- MPK690 OS

With an innovative design that is very easy to put on (like a vest) and with an open design to facilitate diapers changes. Characteristics seamless bonding" technology to protect baby from chafing. Adjustable straps with safety buckles. Allow the abduction of each hip to be regulated independently. Prevent the child from altering them during use. Backstrap. Fits and fastens while maintaining abduction by gravity (not forced).

Ref.	Size
MPK690 OS	One Size



Chest circumference range: 30-55 cm



### Knee support with patella pad and straps - MPK700

Centres and stabilises the kneecap thanks to its removable patella pad. Its spiral side stays stabilise the joint in the mediolateral direction allowing flexion/extension movement. Size 2 features a removable, 100% silicone patella pad which exerts pressure on the patella tendon, relieving pain. Open in the area of the popliteal hollow for greater comfort. Inner fabric with non-slip texture. Fits left and right.

Ref.	Size	Circumference around kneecap (cm)
MPK700 1	T1	20-26
MPK700 2	T2	26-32



### Knee immobiliser - MPK701

Made from a honeycomb technical sports material with soft terry cloth inside to absorb moisture and help maintain healthy skin. Contains two ergonomically shaped side splints and one rear to immobilise the leg in extension. Features an elastic fabric rear panel to provide a snug fit to the back of the knee. Its internal texture prevents slippage. Open kneecap for greater comfort. Fits left and right.

Ref.	Size	Circumference around kneecap (cm)
MPK701 1	T1	24-30
MPK701 2	T2	30-36



### Ankle support - MPK800

Supports and stabilises the joint allowing mediolateral control of the ankle thanks to its front crossover cords. Designed with an open pattern at the back to facilitate fitting and provide comfort when walking. Its velcro fasteners enable compression to be regulated as needed. Fits left and right.

Ref.	Size	Circumference around malleoli (cm)
MPK800 1	T1	16-19
MPK800 2	T2	19-22



# Personalised treatment guide



## Not all patients are the same



Identify the most appropriate orthotic treatment based on the characteristics of your patient and the injury they have.



**A simple practical guide.**  
Flexible and comprehensive help with immediate personalised solutions

With a simple and eminently professional approach, this guide will allow you to apply and understand orthopaedics in a more directed and personal way. Within its pages, you will find all of the necessary information to offer individualised advice to your clients/patients.

**Injury**

**in recovery phase**

**Sport / Work**


**Senior or medium activity**

**Treatment**

**Compression, proprioception, support**

**Compression, proprioception, support**

**Fabric**



**3D elastic**



**2D elastic**

**Line**

**activo sport**

**activo skin**

**Characteristics**

**SkinTech 3D technology**  
• With strategic areas of selective compression  
• AirQuick breathable

**2D flexibility**  
• Elasticity adaptable like a second skin  
• AirQuick breathable


**Which provide**

**Confidence and security for total recovery**

# Personalised Orthotic Treatment Guide : Orthotics

**chronic/  
avoid relapses\***  
Medium activity

Therapeutic heat,  
compression, support,  
stabilisation




Neoprene



- AirLock System technology**
- Keeps warm, expels sweat: breathes
  - Comfortable lightness

**acute**  
High activity

Continuous breathability,  
support,  
stabilisation




Technical-textile fabric



- 100% active breathability**
- Technical-functional pattern
  - Dynamic fitting

**chronic/acute  
(requires immobilisation)**

Immobilisation,  
posture health  
(night-time)



2D technical fabric



- High-frequency technology**
- Seamless manufacture
  - Breathable

**Provide stabilisation and reduce pain without limiting movement**

**Provide immobilisation, reducing pain**

**We move life.**  
We are specialists



\*After major injuries, avoid relapses with Neoprair

### 1.3


#### New line of orthoses for sportspeople

Elastic supports that provide stabilisation and compression to the joints. Preventative during sporting activities. Proprioceptive stimulation. With anatomical shape and seamless for greater comfort without chafing.

#### An effective companion while playing sport




#### AirQuick Fabric



#### Technology that adapts to you

**Aqtivo Sport** elastic supports, which are structured in **lightweight and flexible 3D fibres**, fit like a second skin providing soft and beneficial compression while playing sport. Their strategic selective compression areas protect the muscles and joints from overloads, promote blood circulation and boost proprioception levels for better sports performance.

Designed for maximum breathability, **Aqtivo Sport** elastic supports contain **AirQuick** fabric to quickly absorb and expel moisture to the outside, allowing a continuous flow of air. Dry skin, total comfort.



Display stands :

**activo**  
sport



Counter stand -  
AQTSPXPSI



DISPLAY STAND SIGN

PRIM FLOOR STAND -  
MBSUELOPRIMEXP

Floor stand  
AQTIVO SPORT






**Elastic elbow support - P707**

Lightweight elbow flexure area to facilitate use throughout the day. Anatomical silicone padding to help distribute pressure over the epicondyle. Forearm strap to reinforce pressure at the point of pain. Fits left and right.


Ref.	Size	Elbow circumference (cm)
P707 S	S	23-27
P707 M	M	27-31
P707 L	L	31-35




**Elastic wrist support - P703**

Natural adaptation to the wrist to eliminate strain and relieve pain. Wrist strapping to selectively reinforce compression. Fits left and right.


Ref.	Size
P703	One Size




**Elastic metacarpal wrist support - P704**

Unloads the carpometacarpal joints. Wrist strapping to selectively reinforce compression. Fits left and right.


Ref.	Size	Wrist circumference (cm)
P704 S	S	14-17
P704 M	M	17-20
P704 L	L	20-24




**Elastic thigh support - P709**

Decreasing degree of compression to improve blood circulation, reduce muscle fatigue and promote resistance, minimising the risk of injury. Integrated non-slip system. Fits left and right.


Ref.	Size	Thigh circumference (cm)
P709 S	S	43-48
P709 M	M	48-52
P709 L	L	52-58




**Elastic calf support - P708**

Decreasing degree of compression to improve blood circulation, reduce muscle fatigue and promote resistance, minimising the risk of injury. Integrated non-slip system. Fits left and right.

Ref.	Size	Calf circumference (cm)
P708 S	S	29-35
P708 M	M	35-40
P708 L	L	40-44



## Elastic supports

**aqativo**  
sport



**SKIN  
TECH  
3D**



### Elastic knee support - P700

Lightweight popliteal area to facilitate use throughout the day. Integrated non-slip system. Fits left and right.

Ref.	Size	Knee circumference (cm)
P700 S	S	32-35
P700 M	M	36-39
P700 L	L	40-46
P700 XL	XL	47-52



### Elastic knee support with silicone padding and side stabilisers - P701

Anatomical patella pad to centre the kneecap, helping to dissipate strain. Side stabilising stays for greater mediolateral control. Lightweight popliteal area to facilitate use throughout the day. Integrated non-slip system. Fits left and right.

Ref.	Size	Knee circumference (cm)
P701 S	S	32-35
P701 M	M	36-39
P701 L	L	40-46
P701 XL	XL	47-52



### Elastic ankle support with silicone malleolar pads - P705

Lightweight instep and heel area to facilitate use throughout the day. Anatomical silicone padding over the malleoli to provide stabilisation. Protects against tendon overload and minimises the risk of injury. Fits left and right.

Ref.	Size	Ankle circumference (cm)
P705 S	S	17-20
P705 M	M	20-23
P705 L	L	23-26



### Elastic ankle support with silicone malleolar pads and figure-of-eight strapping - P706

Lightweight instep and heel area to facilitate use throughout the day. Anatomical silicone padding over the malleoli to provide stabilisation. Functional figure-of-eight strapping to enable selective immobilisation and provide stability under high loads. Fits left and right.

Ref.	Size	Ankle circumference (cm)
P706 S	S	17-20
P706 M	M	20-23
P706 L	L	23-26

### 1.4

#### New line of orthoses for everyday use

Elastic supports that provide stabilisation and gentle compression to the joints. Preventative measure during work activities.

Proprioceptive stimulation. With anatomical shape and seamless for greater comfort without chafing.

#### Second skin. Greater protection and well-being

#### AirQuick Fabric



#### Adaptable elasticity. Like a second skin

**Aqtivo Skin** is designed to be a second skin. A line of elastic, discreet and very comfortable supports, manufactured without seams and anatomically shaped to gently protect muscles and joints. Their areas of selective compression and the lightness of their firm but comfortable fibres make **Aqtivo Skin** a range of supports that multiply daily feelings of well-being.

Designed for maximum breathability, **Aqtivo Skin** elastic supports contain **AirQuick fabric** to quickly absorb and expel moisture to the outside, allowing a continuous flow of air. Dry skin, total comfort.



Display stands :

**activo**  
skin



**DISPLAY STAND SIGN**

**PRIM FLOOR STAND -  
MBSUELOPRIMEXP**

Floor stand  
AQTIVO SKIN




Your comfortable second skin



**Elastic wrist support - P703BG**


Natural adaptation to the wrist to eliminate strain and relieve pain. Wrist strapping to selectively reinforce compression. Fits left and right.

Ref.	Size
P703BG	One Size 



**Elastic metacarpal wrist support - P704BG**


Unloads the carpometacarpal joints. Wrist strapping to selectively reinforce compression. Fits left and right.

Ref.	Size	Wrist circumference (cm)
P704BG S	S	14-17 
P704BG M	M	17-20
P704BG L	L	20-24



**Elastic knee support - P700BG**


Lightweight popliteal area to facilitate use throughout the day. Integrated non-slip system. Fits left and right.

Ref.	Size	Knee circumference (cm)
P700BG S	S	32-35 
P700BGM	M	36-39
P700BG L	L	40-46
P700BG XL	XL	47-52



**Elastic knee support with silicone padding and side stabilisers - P701BG**

Anatomical patella pad to centre the kneecap, helping to dissipate strain. Side stabilising stays for greater mediolateral control. Lightweight popliteal area to facilitate use throughout the day. Integrated non-slip system. Fits left and right.

Ref.	Size	Knee circumference (cm)
P701BG S	S	32-35 
P701BG M	M	36-39
P701BG L	L	40-46
P701BG XL	XL	47-52

## Elastic supports


**activo**  
skin



**Elastic elbow support - P707BG**

Lightweight elbow flexure area to facilitate use throughout the day. Anatomical silicone padding to help distribute pressure over the epicondyle. Forearm strap to reinforce pressure at the point of pain. Fits left and right.


Ref.	Size	Elbow circumference (cm)
P707BG S	S	23-27
P707BG M	M	27-31
P707BG L	L	31-35




**Elastic ankle support - P710BG**

Natural adaptation to the ankle to eliminate strain and relieve pain. Lightweight instep and heel area to facilitate use throughout the day. Fits left and right.


Ref.	Size	Ankle circumference (cm)
P710BG S	S	17-20
P710BG M	M	20-23
P710BG L	L	23-26




**Elastic ankle support with silicone malleolar pads - P705BG**

Lightweight instep and heel area to facilitate use throughout the day. Anatomical silicone padding over the malleoli for greater stabilisation. Protects against tendon overload and minimises the risk of injury. Fits left and right.


Ref.	Size	Ankle circumference (cm)
P705BG S	S	17-20
P705BG M	M	20-23
P705BG L	L	23-26




**Elastic ankle support with silicone malleolar pads and figure-of-eight strapping - P706BG**

Lightweight instep and heel area to facilitate use throughout the day. Anatomical silicone padding over the malleoli for greater stabilisation. Functional figure-of-eight strapping to enable selective immobilisation and provide stability under high loads. Fits left and right.

Ref.	Size	Ankle circumference (cm)
P706BG S	S	17-20
P706BG M	M	20-23
P706BG L	L	23-26



1.5

A double effect that  enhances well-being

Our new NEOPRAIR orthoses feature AIRLOCK SYSTEM technology to enable the skin to breathe while maintaining body heat. A dual effect that enhances comfort and pain relief.

**Keep your heat. Breathe**

The fabric, which provides multidirectional elasticity and three-dimensional compression, makes NEOPRAIR orthoses 100% lightweight and adaptable to the shape and mobility of the joints, enhancing its ability to offer support and stabilisation during sports activities.



**Expels sweat**

**Breathability**

Its AirLock System technology generates channels in the neoprene which enable sweat molecules to evaporate allowing the skin to breathe.

**AirLock SYSTEM** 



**Keep your heat**

**Therapeutic heat**

The internal heat cycle is maintained thanks to the intrinsic properties of neoprene. Heat that relieves pain and reduces recovery time.

## Neoprene supports

**NEOPRENE  
SUPPORTS**  
100% lightweight

breathable  
neoprene  
**AirLock SYSTEM**



**DISPLAY STAND SIGN**

**PRIM FLOOR STAND -  
MBSUELOPRIMEXP**

Floor stand  
NEOPRAIR / MOTTIO /  
AIRMED / MYPRIM  
KIDS / AIRTEX





**Knee support - NPOS136**

Patella pad to stabilise the kneecap. 4 splints to improve lateral stability. Inner silicone bands with anti-slip effect. Opening in the popliteal area for greater comfort. Double fastening strap for an optimum fit. Bilateral.

Ref.	Size
NPOS136	One size



**Elbow support- NPOS173**

Double fastening strap for an optimum fit. Enables control of pressure. Bilateral.

Ref.	Size
NPOS173	One size



**Ankle support - NPOS111**

Combines breathable neoprene of different thicknesses to provide greater comfort in the foot and fastening area. Bilateral.

Ref.	Size
NPOS111	One size



**Thigh support - NPOS147**

Inner silicone bands with anti-slip effect. Double fastening strap for an optimum fit. Bilateral.

Ref.	Size
NPOS147	One size

ONE  
SIZE **S**

Neoprene supports



neoprair



Calf support - NPOS120

Combines breathable neoprene of different thicknesses to provide greater comfort in the calf. Inner silicone bands with anti-slip effect. Triple fastening strap for an optimum fit. Bilateral.

Ref.	Size
NPOS120	One size



Patellar strap - NPOS132

Silicone patellar pad that unloads the patellar tendon. With non-slip effect. Bilateral.

Ref.	Size
NPOS132	One size



Epicondylitis strap - NPOS172

Silicone condyle pad that presses on the tendon to help eliminate pain. Inner silicone band to prevent the orthosis from rotating. Bilateral.

Ref.	Size
NPOS172	One size



Bandage-style wrist support - NPOS164

Easily adjustable non-elastic strap. Bilateral.

Ref.	Size
NPOS164	One size



**Wraparound wrist support - NPOS167**

Wraps around the trapeziometacarpal and wrist joint. Bilateral.

Ref.	Size
NPOS167	One size



**Wraparound wrist and thumb support - NPOS165**

Wraps around the trapeziometacarpal and wrist joint. Reduces pain and instability of the thumb. Bilateral.

Ref.	Size
NPOS165	One size



**Lumbosacral band - NPOS154**

Large fastening area that makes fitting easier and improves comfort.

Front height: 24cm.  
Back height: 24 cm.

Ref.	Size
NPOS154	One size



**Lumbosacral back support - NPOS155**

4 stays that provide greater stability in the lumbar region. Anatomical fastening for an improved fit.

Front height: 16 cm.  
Back height 24 cm.

Ref.	Size
NPOS155	One size

## Display Stand

mottio<sup>HF</sup>  
sport range

1.6



### Keep moving

Designed to alleviate impacts to the lateral/medial epicondyles and patellar tendon during sporting activity. Its anatomical design and soft-to-the-touch fabrics adapt to the user's movement, enabling him to perform his activities with total freedom to achieve maximum performance.

*Breathable seamless fabric for maximum comfort.*

### Counter stand - EXPMOTTIO

Measurements  
35x48 cm.



### Epicondylitis strap - T29

With 100% silicone anatomical condyle pad with massage effect and repositionable. Elastic strap for dynamic compression adjustment. Anatomical structure to make adjustment to the forearm easier. Hypoallergenic non-slip internal support strap. Latex-free.

Ref.	Size	Side
T29	One size	Fits left and right



Anatomical condyle pad

## Supports straps for high performance



Anatomical patellar pad



### Infrapatellar strap - T129

With 100% silicone anatomical patellar pad with massage effect. Non-slip. Moulded pad in the area of the popliteal hollow to provide greater comfort. Double elastic strap for dynamic compression adjustment. Latex-free.

Ref.	Size	Side
T129	One size	Fits left and right



### Breathable textile technology for healthy skin

Fabric supports made from 100% breathable state-of-the-art technical fabric which provides high elasticity to the item and heat, compression and stabilisation to the treated area, keeping the skin healthy thanks to its high breathability.

**Get moving and let your skin breathe.**



Multi-cell textile design for greater comfort.



Highly-elastic fabric capable of adapting to dynamic changes.




Anatomical seams for a comfortable fit.



One, two or three fastening straps for easy adjustment and rotation prevention.




One size design that adapts to any contour.



Comfort zone

Highly-elastic and breathable technical impact fabric that adapts to the dynamics of movement, providing the patient with total comfort.



Compression zone

Breathable, elastic, wraparound velour fabric that enables the compression required by the patient to be adjusted at any time.



DISPLAY STAND SIGN

PRIM FLOOR STAND -  
MBSUELOPRIMEXP

Floor stand  
NEOPRAIR / MOTTIO /  
AIRMED / MYPRIM  
KIDS / AIRTEX



Get moving and let your skin breathe





**Elbow support - OST273**

Anatomical design for a better fit to the joint area. With open pores to maximise ventilation of the elbow.  
Fits left and right.

Ref.	Size
------	------

OST273 One size



**Epicondylitis strap - OST272**

Foam pad that exerts pressure on the epicondyle to relieve pain.  
Wraparound fastening for easy individualised adjustment.  
Fits left and right.

Ref.	Size
------	------

OST272 One size



**Wraparound wrist support - OST265**

With wraparound fastening that allows for quick and easy fitting.  
Fits left and right.

Ref.	Size
------	------

OST265 One size

## Fabric supports

AirTox  
Breathable

ONE  
SIZE  
TALLA  
ÚNICA



### Wrist and thumb support - OST264

Flexible design that stabilises the thumb and enables the functionality of the hand to be maintained, relieving pain. Left and right-hand versions.

Ref.	Size	Side
OST264	One size	L-Left
OST264		R-Right



### Knee support - OST218

With triangular anatomical patellar opening. Comfortable interior foam pad. Open popliteal area for greater freedom of movement. Medial-lateral stabilisers. Fits left and right.

Ref.	Size
OST218	One size



### Patellar strap - OST232

Silicone pad that exerts pressure on the patellar tendon to relieve pain. Fits left and right.

Ref.	Size
OST232	One size



**Thigh support - OST247**

Featuring an anatomical design that provides high compression and support. Interior silicone strip to prevent slippage. Fits left and right.

Ref.	Size
OST247	One size



**Calf support - OST220**

Anatomical design that ensures a snug fit to the calf muscle, providing greater compression and support. Fits left and right

Ref.	Size
OST220	One size



**Cross-strap ankle support - OST212**

Figure-of-eight strapping for greater and more effective bracing. Open heel to enable joint movement. Fits left and right.

Ref.	Size
OST212	One size

## Fabric supports

AirTex  
Breathable

ONE  
SIZE  
TALLA  
ÚNICA




### AIRTEX lumbar support belt - OST255BG

Intelligent modular design that adapts to any size. Lumbar support with 4 removable and adjustable steel stays. Adjustable and removable elastic tensioning bands with dual fastening system. With Velcro tabs to facilitate attachment and detachment of each piece. Moderate compression. Breathable. Latex-free.

Front height: 16 cm Back height: 26 cm

Ref.	Size	Colour	Waist circumference (cm)
OST255BG	One size	Beige	80-125




### AIRTEX lumbar support belt - OST255G

Intelligent modular design that adapts to any size. Lumbar support with 4 removable and adjustable steel stays. Adjustable and removable elastic tensioning bands with dual fastening system. With Velcro tabs to facilitate attachment and detachment of each piece. Moderate compression. Breathable. Latex-free.

Front height: 16 cm Back height: 26 cm

Ref.	Size	Colour	Waist circumference (cm)
OST255G	One size	Grey	80-125





**Latest generation wristbands  
for the most demanding activity.**

With perforated ventilation areas to offer maximum breathability.

Its open and seamless design prevents friction and irritation maximizing comfort.

tecnología de  
alta frecuencia  
high frequency technology



**100% LATEX FREE** ✓



## Immobilizers wrist support

airMed®



### Thumb abduction orthosis – AM202G

Immobilises the trapeziometacarpal and metacarpophalangeal area of the thumb in a functional position thanks to its anatomical splints made from malleable aluminium. With velcro fastening for an optimum fit. Fits left and right.

Ref.	Size	Wrist circumference (cm)	Thumb circumference (cm)
AM202G S	S	14-18	5-7
AM202G L	L	18-22	6-8

Also available in beige - AM202



### Immobilising wrist brace with palmar splint - AM205G

Features a palmar plus splint made from malleable aluminium which allows dorsal and palmar flexion to be limited, helping to reduce inflammation of the median nerve. With 2 aluminium dorsal stays that limit the lateral movements of the wrist. Features an adjustable closure band for extra individualised compression of the wrist. Fits left and right.

Ref.	Size	Side	Length (cm)	Wrist circumference (cm)
AM205G S	S	Fits left and right	20.5	14-18
AM205G L	L	Fits left and right		18-22



### 1.8



**Short immobilising wrist brace - AM203G**

Limits dorsal and palmar flexion, helping to reduce inflammation of the median nerve thanks to its malleable aluminium palmar splint. Its 4 dorsal stays limit the lateral movements of the wrist. Features an adjustable closure band for extra individualised compression of the wrist.  
Left and right versions

Ref.	Size	Side	Length (cm)	Wrist Circumference (cm)
AM203G SR	S	R-Right	15	14-18
AM203G SL	S	L-Left		14-18
AM203G LR	L	R-Right		18-22
AM203G LL	L	L-Left		18-22



Also available in beige - AM203



**Long immobilising wrist brace - AM204G**

Ref.	Size	Side	Length (cm)	Wrist circumference (cm)
AM204G SR	S	R-Right	19.5	14-18
AM204G SL	S	L-Left		14-18
AM204G LR	L	R-Right		18-22
AM204G LL	L	L-Left		18-22



Also available in beige - AM204



**Thumb immobilising attachment – AM201G**

Its malleable aluminium splint immobilises the thumb area, relieving the pressure of the synovial sheath on the short and long abductor tendons. For use with items AM203G, AM204G and AM205G.  
Fits left and right.

Ref.	Size
AM201G	One size

Also available in beige - AM201

# Immobilizers wrist support




**Short wrist Support With thumb  
- AM206G**




**Long wrist Support With thumb  
- AM208G**

Breathable, lightweight and adaptable fabrics With perforated ventilation areas. Functional immobilization splints with anatomically contoured aluminium splints. "itself" fit system Adjust the orthosis easily and safely using just one hand

Ref.	Size	Side	Length (cm)	Wrist Circumference (cm)
AM206G SR	S	R-Right	16	14-18
AM206G SL	S	L-Left		14-18
AM206G LR	L	R-Right		18-22
AM206G LL	L	L-Left		18-22



Ref.	Size	Side	Length (cm)	Wrist circumference (cm)
AM208G SR	S	R-Right	20	14-18
AM208G SL	S	L-Left		14-18
AM208G LR	L	R-Right		18-22
AM208G LL	L	L-Left		18-22





1.9

**An indispensable range for your shelves**

The ESSENCIAL line of stabilising wrist supports features mouldable and removable aluminium splints in the palmar and/or thumb area for individualised, comfortable and effective stabilisation of the joints.

An essential item that you need to have on your shelves.



**Stabilising wrist support with V-shaped splint - C400**

With 2 removable and mouldable aluminium splints positioned in the shape of a V for greater comfort. With secure Velcro fasteners and easy opening. Internal micro-massage effect fabric. Breathable. Latex-free.

Fits left and right.

Ref.	Size	Length (cm)	Wrist circumference (cm)
C400 S	S	6	12-15
C400 M	M	6	15-17
C400 L	L	6	17-19
C400 XL	XL	6	19-21



**Stabilising wrist support with thumb splint - C700**

With removable and mouldable aluminium thumb splint. With secure Velcro fasteners and easy opening. Internal micro-massage effect fabric. Breathable. Latex-free.

Ref.	Size	Side	Length (cm)	Wrist circumference (cm)
C700 SR	S	R-Right	13	12-15
C700 SL	S	L-Left	13	12-15
C700 MR	M	R-Right	13	15-17
C700 ML	M	L-Left	13	15-17
C700 LR	L	R-Right	13	17-19
C700 LL	L	L-Left	13	17-19
C700 XLR	XL	R-Right	13	19-21
C700 XLL	XL	L-Left	13	19-21



# Stabilising elastic wrist support

ESSENCIAL



## Stabilising wrist support with thumb and palmar splints - C800

With removable and mouldable aluminium splints in the thumb and palmar area. With secure Velcro fasteners and easy opening. Internal micro-massage effect fabric. Breathable. Latex-free.

Ref.	Size	Side	Length (cm)	Wrist circumference (cm)
C800 SR	S	R-Right	20	12-15
C800 SL	S	L-Left	20	12-15
C800 MR	M	R-Right	20	15-17
C800 ML	M	L-Left	20	15-17
C800 LR	L	R-Right	20	17-19
C800 LL	L	L-Left	20	17-19
C800 XLR	XL	R-Right	20	19-21
C800 XLL	XL	L-Left	20	19-21



## Short stabilising wrist support with palmar splint - C500

With removable and mouldable aluminium palmar plate. With secure Velcro fasteners and easy opening. Internal micro-massage effect fabric. Breathable. Latex-free.

Ref.	Size	Side	Length (cm)	Wrist circumference (cm)
C500 SR	S	R-Right	18	12-15
C500 SL	S	L-Left	18	12-15
C500 MR	M	R-Right	18	15-17
C500 ML	M	L-Left	18	15-17
C500 LR	L	R-Right	18	17-19
C500 LL	L	L-Left	18	17-19
C500 XLR	XL	R-Right	18	19-21
C500 XLL	XL	L-Left	18	19-21



## Long stabilising wrist support with palmar splint - C600

With removable and mouldable aluminium palmar plate. With wrist strap to reinforce the support. With secure Velcro fasteners and easy opening. Internal micro-massage effect fabric. Breathable. Latex-free.

Ref.	Size	Side	Length (cm)	Wrist circumference (cm)
C600 SR	S	R-Right	22	12-15
C600 SL	S	L-Left	22	12-15
C600 MR	M	R-Right	22	15-17
C600 ML	M	L-Left	22	15-17
C600 LR	L	R-Right	22	17-19
C600 LL	L	L-Left	22	17-19
C600 XLR	XL	R-Right	22	19-21
C600 XLL	XL	L-Left	22	19-21





**Back support with 2 bands**  
- PRS602BG

Band made of breathable stretch fabric, soft and pleasant to the touch. Dual straps that adapt to the different areas of the abdomen, providing better comfort and compression. Its wide micro velcro fastening area allows for a secure and individualised fit. It uses "high frequency technology" for a 100% seamless design, more comfortable and better for the skin. It is designed both for abdominal and costal injuries. Latex free.

Size	Circumference (cm)	Height (cm)
S	70-90	16
M	90-110	16
L	110-130	16
XL	130-145	16



**Back support with 3 bands**  
- PRS603BG

Band made of breathable stretch fabric, soft and pleasant to the touch. Three straps that adapt to the different areas of the abdomen, providing better comfort and compression. Its wide micro velcro fastening area allows for a secure and individualised fit. It uses "high frequency technology" for a 100% seamless design, more comfortable and better for the skin. Latex free.

Size	Circumference (cm)	Height (cm)
S	70-90	24
M	90-110	24
L	110-130	24
XL	130-145	24



**Back support with 4 bands**  
- PRS604BG

Band made of breathable stretch fabric, soft and pleasant to the touch. Four straps that adapt to the different areas of the abdomen, providing better comfort and compression. Its wide micro velcro fastening area allows for a secure and individualised fit. It uses "high frequency technology" for a 100% seamless design, more comfortable and better for the skin. Latex free.

Size	Circumference (cm)	Height (cm)
S	70-90	32
M	90-110	32
L	110-130	32
XL	130-145	32

# Back Supports



## MOVE Cotton stockinette back support- PRS600

Closed stockinette support made from knitted fabric. Adapts perfectly to all anatomies, achieving maximum discretion. Offers optimum freedom of movement. Latex free.

## MOVE Cotton stockinette with open - PRS601

Lumbosacral back support made from knitted fabric. Provides light support and allows maximum freedom of movement. Adapts perfectly to all anatomies, achieving maximum discretion. Latex free.

## MOVE+ Back support - PRS606

Lumbosacral back support made from breathable monofilament to keep the skin dry and free from moisture caused by perspiration. With two plastic stays to provide stability in the lumbar area. Latex free.

## MOVE+ Reinforced back support - PRS607

Lumbosacral back support made of dense elastic fabric. With four plastic stays for enhanced support in the lumbar area. Anatomically shaped design with a narrower front for improved comfort. Latex free.



Size	Circumference (cm)	Height (cm)
XS	65-75	26
S	75-95	29
M	95-115	29
L	115-130	29



Size	Circumference cintura (cm)	Height (cm)
S	75-95	26
M	95-115	29
L	115-130	29



Size	Circumference cintura (cm)	Height (cm)
S	75-95	24
M	95-115	24
L	115-130	24
XL	130-140	24



Size	Circumference cintura(cm)	Height frontal(cm)	Height espalda(cm)
S	70-90	14	24
M	90-120	14	24
L	120-140	14	24




### Support levels

The elements and materials we use in the design of each orthosis: fabrics, tensioners, stays, thermoplastics, metals, etc., as well as the shape, size and thickness we apply, allows us to customise the treatment by establishing different support levels.



Factors involved in good spine support:

- Lumbar compression
- Lumbar support
- Abdominal compression
- Shoulder retropulsion (at dorsal level)



Exclusive double-channel boostfit fabric

**BOOST fit**

- + Fitting
- + Compression
- + Comfort



# Lumbosacral back supports



**MAMMY Maternity back support - PRS615**

Front anatomical band made of breathable, soft and light velour, which raises the abdomen, relieving the feeling of heaviness in the pelvic area. It incorporates "boostFit" breathable stretch fabric, with a double channel structure that improves adaptation to the body, providing better compression and comfort. 2 lumbar containment stays that support the spine and improve posture, relieving back pain. 2 lateral tensioners that can be adjusted to unload pressure depending on the needs of the mother-to-be. It uses seamless, "high-frequency technology", which promotes discretion and comfort. With help guides for the correct fitting of the belt in each pregnancy stage. Latex free.




Color	Size	Height front (cm)	Height back (cm)
Beige	One Size	7	13
Grey	One Size	7	13



**Semi-rigid lumbosacral back support - PRS625N**

Containment belt for work, made of breathable elastic lycra tulle, with a velcro front fastening. Rear support consisting of 4 stays that provide greater stability in the lumbar region. Plus 2 inguinal reinforcement stays. Duo Control system with a double tensioning band lets you adjust the lumbar support in line with your work-related needs. It has non-slip inner bands. It can be fitted with shoulder straps for greater safety (supplied separately as an accessory Reference: PRS6ATR OS).


Size	Circumference (cm)	Height front (cm)	Height back (cm)
XS	75-85	15	24
S	85-95	15	24
M	95-105	15	24
L	105-115	15	24
XL	115-125	15	24
XXL	125-135	15	24
XXXL	135-150	15	24



**LABORAL+ Semi-rigid lumbosacral back support- PRS635N**

Containment belt for work, made of breathable monofilament stretch fabric offering great compression, with a Velcro front fastening. Rear support consisting of 4 steel stays that provide greater containment in the lumbar region. Reinforced with 2 inguinal stays. Duo Control system with a double tensioning band lets you adjust the lumbar support in line with your work-related needs. Handles on the fastening and on each tensioner for easy fitting. It incorporates internal padded back straps to maximise comfort.

Size	Circumference (cm)	Height front (cm)	Height back (cm)
S	85-95	17	26
M	95-105	17	26
L	105-115	17	26
XL	115-125	17	26
XXL	125-135	17	26





Support level

## CARE Semi-rigid lumbosacral back support - PRS610

Semi-rigid back support belt made of “boostFit” stretch fabric, highly breathable, with a double channel structure that enhances the adaptation to the body, providing better compression and comfort. It has a rear crossover for greater support and 4 thick stays, which provide great containment in the lumbar area. The lumbar crossover is the most suitable system when a high degree of support is required in a specific area. Micro velcro fastening for a secure fit, and front finger loops to facilitate adaptation in cases of osteoarthritis in the hands or other disorders. It uses seamless, “high-frequency technology”, which promotes discretion and comfort. Latex free.



Exclusive double-bannel fabric

Size	Circumference (cm)	Height front (cm)	Height back (cm)
XS	70-80	17	26
S	77-90	17	26
M	87-100	17	26
L	95-110	17	26
XL	107-120	17	26
XXL	117-130	17	26



Also available in beige



Support level

## CARE+ Semirigid lumbosacral back support - PRS620

Semi-rigid back support belt made of “boostFit” stretch fabric, highly breathable, with a double channel structure that enhances the adaptation to the body, providing better compression and comfort. It has 4 thick stays and double criss-cross tensioner in the lumbar area with front fastening, providing stability and better containment. It also has a foam lumbar plate that provides heat and comfort to the area and can be removed when not required. Micro velcro fastening for a secure fit, and front finger loops to facilitate adaptation in cases of osteoarthritis in the hands or other disorders. It uses seamless, “high-frequency technology”, which promotes discretion and comfort. Latex free.



Exclusive double-bannel fabric

Size	Circumference pelvis (cm)	Height front (cm)	Height espalda (cm)
XS	70-80	17	26
S	77-90	17	26
M	87-100	17	26
L	95-110	17	26
XL	107-120	17	26
XXL	117-130	17	26



Also available in beige







**EXCELLENT sling - 909**

With rubber shoulder support. With thumb loop. Height adjustable with buckles. Made from high-quality fabric. Breathable. Latex-free. Fits left and right.

Ref.	Size	Forearm and palm of hand length (cm)
909 S	S	< 30
909 M	M	30-42
909 L	L	> 42



**PRIM forearm sling - 901**

Armsling made in breathable fabric offering high comfort. Designed with a regulation system by buckle and Velcro® which allows adjustment to the most convenient position to the patient. Breathable. Latex-free. Fits left and right.

Ref.	Size	Forearm and palm of hand length (cm)
901 S	S	28 - 36
901 M	M	35 - 42
901 L	L	41 - 49
901 XL	XL	48 - 56



**PRIM children's sling - 907**

Easy placement children's sling. Made from polyester and cotton. Fits left and right.

Ref.	Size	Measurements
907	One size	28x14 cm

## Others



### Armsling/ shoulder immobilizer - 903

Made in a padded fabric. This armsling is adjustable in height by means of buckle and Velcro®. Designed with an immobilization band.


Ref.	Size	Length of forearm and palm (cm)
903 S	S	28-36
903 M	M	35-42
903 L	L	41-49
903 XL	XL	48-56



### CONFORT shoulder support - 335

Promotes retropulsion of the shoulders and keeps the back straight, helping to correct kyphotic postures. With a paravertebral band for a controlled fit. With removable anatomical pads in the armpit area. Breathable elastic monofilament fabric. Latex-free.


Ref.	Size	Circumference below the chest (cm)
335 S	S	75-93
335 M	M	88-107
335 L	L	102-120




### ELCROSS LUX shoulder support - 3138BG

Orthosis made of breathable double layer lycra material. It uses carbon fibre with Resistex® Carbon technology, especially developed to offer to the patient maximum comfort and breathability. Maintains the spine in anatomical position, avoiding incorrect and/or kyphotic positions. It has two rear stays that help to improve postural correction. Its straps promote retropulsion of the shoulders and are reinforced with crossed elastic straps in the dorsal area, ensuring adequate correction. Latex-free.

Ref.	Size	Circumference below the chest (cm)
3138BG S	S	75-93
3138BG M	M	88-107
3138BG L	L	102-120
3138BG XL	XL	
3138BG XXL	XXL	




Antistatic  
Thermoregulator  
Protective (UV, static energy)

# Cervical Collar



## Straight soft cervical collar - CC19

Collar made of foam covered by cotton stockinette. Breathable. Straight collar. Rear reinforcement.

Ref.	Size	Measurement (cm)*	Height (cm)
CC19 XS	XS/low	32 - 44	8
CC19 S	S/low	37 - 49	8
CC19 M	M/high	42 - 54	10
CC19 XL	L/high	48 - 60	10

\* Circumference and height of the neck




## Contoured soft cervical collar - CC20

Collar made of foam covered by cotton stockinette. Breathable. Contoured collar for a better fitting and comfort. Rear reinforcement. Beige colour.

Ref.	Size	Circumference (cm)*	Height (cm)
CC20 PB	S/low	42 - 58	8
CC20 PM	S/medium	42 - 58	9
CC20 PA	S/high	42 - 58	10
CC20 GB	L/low	47 - 63	8
CC20 GM	L/medium	47 - 63	9
CC20 GA	L/high	47 - 63	10

\* Circumference and height of the neck



## Semirigid contoured cervical collar - CC121

Collar made of polyurethane foam and covered by cotton stockinette. Contoured collar with a plastic reinforcement of 1 mm polyethylene, providing more rigidity. This reinforcement can be removed depending on the rigidity the patient may need. Beige colour.

Ref.	Size	Circumference (cm)*	Height (cm)
CC121 PB	S/low	35 - 38	8
CC121 PM	S/medium	35 - 38	9
CC121 MB	M/low	38 - 41	8
CC121 MM	M/medium	38 - 41	9
CC121 GM	L/medium	41 - 43	9
CC121 GA	L/high	41 - 43	11
CC121 XLM	XL/medium	43 - 48	9
CC121 XLA	XL/high	43 - 48	11

\* Circumference and height of the neck



## Post-surgical shoes


### PRIM's post-surgical shoes

feature a Velcro closure system that can adjust to feet with significant deformities or post-operative bandages or plaster casts.

Indicated for post-surgery of the foot and toes, plaster casts or immobilising bandages and swelling that prevents the use of conventional shoes.



non-slip sole



fits left and right



extensive size range



### Open post-surgical shoe - PS100

Easy and secure fastening system. Non-slip sole. Made from breathable fabric to ensure comfortable use. Fits left and right.


Ref.	Size	Shoe size
PS100 0	(0) XS	33/35
PS100 1	(1) S	36/38
PS100 2	(2) M	39/40
PS100 3	(3) L	41/42
PS100 4	(4) XL	43/45
PS100 5	(5) XXL	46/48



### Open post-surgical talipes calcaneus shoe - PS200

Easy and secure fastening system. With extra ankle support. Non-slip sole. Made from breathable fabric to ensure comfortable use. Fits left and right.

Ref.	Size	Shoe size
PS200 0	(0) XS	33/35
PS200 1	(1) S	36/38
PS200 2	(2) M	39/40
PS200 3	(3) L	41/42
PS200 4	(4) XL	43/45



A close-up photograph of a human foot, specifically the big toe area. A white, crescent-shaped protective strip is applied to the side of the big toe. The foot is resting on a light-colored wooden surface. The background is softly blurred, showing the rest of the foot and the lower leg.

# 2

- 2.1. Comforsil
- 2.2. Comforgel
- 2.3. Plantillas Comforgel
- 2.4. OrthoPrim

## 2.1

### Silicone podiatry products

A new range of Comforsil silicone podiatry products that offer a practical and high-quality solution to foot problems.

Products made from 100% pure, hypoallergenic and non-toxic silicone.

Innovative, practical and quality solutions for problems associated with the feet.

Strength. The properties of silicone ensure long durability and crack resistance regardless of how long and intensively the products need to be used.

Comfort. They gently adapt to the feet ensuring recovery and cushioning at loading points and in pressure areas.



**Floor stand -  
CS/CG EXP 4G-2020**

Measurements  
47 x 28 x 20 cm



**Floor stand -  
CS S/C EXP 4G-2020**

Measurements  
47 x 28 x 20 cm





Toes



**Big toe separator (hourglass-shaped) - CC218**

**Extra-thin toe separator (crescent-shaped) - CC206**

**Toe crests with loop - CC226**

**Protective toe cap - CC227**

**Portable manicure set - CC284**

Straightens and aligns the big toe. Eases bunion pain.

Prevents friction between the toes. Placement on any toe.

Correct toe flexion.

Toe and nail protection  
Reinforcement at the tip of the toe.

Nail care and beauty.  
Includes a set of 5 drills (for polishing, filing and shining the nails).  
Works with 2 batteries (not included).  
Do not use rechargeable batteries.

x2

Ref.	Sizes
CC218S	S
CC218M	M
CC218L	L

x3

Ref.	Sizes
CC206S	S
CC206M	M
CC206L	L

x2

Ref.	Sizes
CC226S	S
CC226M	M
CC226L	L

x2 Anti-chafing

Ref.	Sizes
CC227S	S
CC227L	L

Ref.  
CC284

x1 + SET





Toes



**Tubular toe bandage - PRO0108**

Soft toe protection against chafing, calluses and blisters. Cotton-lined interior.



Ref.	Sizes
PRO0108S	S
PRO0108M	M
PRO0108L	L



**Night-time bunion corrector - CC250/1/2**

Correctly positions the big toe. Eases pain during the night. Avoid walking with the bunion corrector in place.

Ref.	Sizes
CC250	One size
CC251	Right (x1)
CC252	Left (x1)



**Elastic fabric protector - CC264**

Elastic protector with interior silicone pad. Made from soft, seamless elastic material with different compression zones.



Ref.	Sizes
CC264S	S
CC264L	L



**Bunion corrector for daily use - CC450**

Daytime hallux valgus corrector with protective silicone pad in the bunion area. Made from soft, seamless elastic material with different compression zones.



Ref.	Side
CC450R	R
CC450L	L



Toes



Repositioning toes - CC510

Tapes to reposition the fingers of the feet in case of fractures or overlaps. Light and comfortable, it can be used with any type of footwear. Washable and reusable.

Ref.	Sizes	Width Long ring overall length
CC510M	M	2x4x9cm
CC510L	L	3x5x11,5cm



Forefoot



Metatarsal pads with 2nd-toe holding loop - CC220

Shock-absorbing (high-heeled shoes). Anti-slip recesses.

Ref.	Sizes
CC220	One size



Plantar protector with pad - CC256

Plantar protector made from soft, seamless elastic fabric. Contains non-slip material in the instep area to prevent movement during use. Features an extra-soft silicone pad inside.

Ref.	Sizes
CC256S	S
CC256L	L



Padded elastic protector - CC320

Triple action elastic protector with interior 100% silicone pad. Made from seamless elastic fabric with different compression zones. Relieves footwear pressure on bunions, metatarsals and calluses.

Ref.	Sizes
CC320SP	SP
CC320LP	LP



SEAMLESS

Forefoot protectors with pad - CC225

Contoured forefoot protector made from seamless elastic fabric. Contains non-slip material in the instep area to prevent movement during use. Includes an extra soft silicone pad inside that covers the area from the tips of the toes to the metatarsals.

Ref.	Sizes
CC225S	S ≥ 22 cm
CC225L	L 22 - 26 cm



Post-surgery elastic metatarsal support - CC322

Adequate forefoot support. Perfect fit and support.

Ref.	Sizes
CC322S	S
CC322M	M
CC322L	L



Elastic splay foot bandage with pad - CC229

Elastic, maximum support and fit. Pad for sole and metatarsal relief.

Ref.	Sizes
CC229S	S
CC229M	M
CC229L	L



Mini insoles for rest - CC257

Comfort and cushioning. Anti-slip recesses.

Ref.	Size
CC257	One size



Forefoot



Lined mini insoles for rest - CCF401

Comfort and cushioning.  
Anti-slip recesses.  
Soft and comfortable lining.



Ref.	Size
CCF401	One size



Heel



Self-adhesive heel protector - CC213

Heel protection.  
Prevents rubbing.



Ref.	Size
CC213	One size



Calcaneal spur heel pads with removable plug - CC215

Pronation corrector.  
Unloads the plantar tendon insertion.



Ref.	Sizes
CC215S	S (34-37)
CC215M	M (38-40)
CC215L	L (41-43)



Heel cups - CC212

Less thickness - better fit .  
Lower density silicone in the heel.



Ref.	Sizes
CC212S	S (34-37)
CC212M	M (38-40)
CC212L	L (41-43)




Heel



**Non-slip heel pads for rest - CC214**

Less thickness - greater fit to footwear.  
Anti-slip recesses.



x2  

Ref.	Sizes
CC214S	S (34-37)
CC214M	M (38-40)
CC214L	L (41-43)



**Lined heel pads for rest - CCF214**

Less thickness - better fit.  
Anti-slip recesses.

x2  

Ref.	Sizes
CCF214S	S (34-37)
CCF214M	M (38-40)
CCF214L	L (41-43)



**Silicone heel lifts - CC247**

Reduce loads on the ankle, hip and spine.

x2  

Ref.	Sizes
CC247S	S (3mm)
CC247M	M (5mm)
CC247L	L (8mm)



Sole



Extra-thin lined insoles - CCF313

Relief in the metatarsal area. Cushion impacts on the heel.



Ref.	Sizes
CCF313S	S (35-36)
CCF313M	M (37-38)
CCF313L	L (39-40)
CCF313XL	XL (41-42)



Insoles for tired feet and legs - CC203

Areas of softer silicone, greater comfort. Relief on ligaments and muscles.



Ref.	Sizes
CC203S	S (35-39)
CC203M	M (40-43)
CC203L	L (44-46)



Insoles with shock-absorbing cushioning - CC209

Maximum reinforcement for the heel and metatarsals. Areas of softer silicone, greater comfort.



Ref.	Sizes
CC209S	S (35-39)
CC209M	M (40-43)
CC209L	L (44-46)



Sole



'Odour Control' charcoal insoles - CC297

Absorb and neutralise bad odour.  
3 technical layers.



Ref.	Size
CC297	One size



Plasters



'Aqua' self-adhesive bunion patches - CC330

Distribute pressure and reduce friction.



Ref.	Sizes
CC330	One size



'Aqua' adhesive callus patches - CC327

Distribute pressure and reduce friction.



Ref.	Sizes
CC327	One size



Plasters



'Aqua' adhesive trimmable strip - CC329

Distribute pressure and reduce friction.



Ref.	Sizes
CC329	9,5x9,5cm



Adhesive trimmable felt strip - CC249

Protection against irritation. Ideal for the area between the toes.



Ref.	Sizes
CC249	22 x 7cm



'Aqua' shin protector - CC332

Distributes pressure. Decreases rubbing.



Ref.	Sizes
CC332	9,5x9,5cm



Expositor mostrador  
- CS/CG EXP 4G-2020

Medidas  
47 x 28 x 20 cm

### Gel podiatry products

ComforGel is a new range of podiatry products which offers comfortable solutions perfectly adaptable to all areas of the foot that need treatment.

Products made from 100% pure gel.

Designed to provide maximum comfort and flexibility by protecting the feet where special care is required.

High flexibility and maximum adaptability. They can adapt to the areas to be treated ensuring greater absorption of impacts, reducing friction and relieving discomfort caused in the feet.





Toes



### Protective tubes with inner gel lining

Protect injured or irritated areas of the toes. Relieve the pressure of footwear and protect against friction and chafing



Ref.	Measurements	Details
CC245	Ø1,5cm/0,59" x 2	Mesh
CC325	Ø1,5cm/0,59" x 2	Canale
CC253	Ø2,5cm/0,98" x 2	Canale
PRO5016	Ø1,5cm/0,59" x12	Mesh
PRO5018	Ø1,5cm/0,59" x12	Canale
PRO5017	Ø2,5cm/0,98" x12	Canale



### Toe cap with elastic fabric - CC267

Made from gel and lined with elastic fabric. Protects the toes from impacts and chafing.

Ref.	Sizes
CC267S	S
CC267L	L



### Toe protectors covered with elastic fabric - CG1201

Toe protectors with gel disc that alleviate discomfort caused by calluses, hammer toes and toe injuries.

Ref.	Sizes
CG1201SM	S/M
CG1201LXL	L/XL



### Protective toe cap - CG0901

Protects the toes from impacts and chafing.

Ref.	Sizes
CG0901SM	S/M
CG0901L/XL	L/XL



### Bunion protector with hourglass-shaped separator - CG2301

Relieves bunion pain and helps preserve joint alignment.

Ref.	Sizes
CG2301	One size



**Thick toe separator (hourglass-shaped) - CG1101**

x2

Extra soft to straighten and align the big toe, protecting against friction and overlapping toes.

Ref.	Sizes
CG1101S	S
CG1101M	M
CG1101L	L



**Bunion protector with adjustable loop - CG1801**

x1

Soft and flexible protector that protects bunions from the pressure and friction of footwear.

Ref.	Sizes
CG1801	One size



**Toe separator with loop - CG2101**

x2

Used for the protection of corns and prevention of interdigital rubbing, it helps toe alignment and reduces friction.

Ref.	Sizes
CG2101M	M
CG2101L	L



**Toe crest with support loop - CG0801**

x1

Relief for patients with claw or hammer toes.

Ref.	Sizes
CG0801	One size



**Little toe bunion protector - CG1901**

x1

Protects bunions on the little toe from the pressure and friction of footwear.

Ref.	Sizes
CG1901	One size



**Elastic fabric bunion protector - CG1601**

x1

Relieves bunion pain caused by the pressure of footwear and reduces the risk of discomfort from rubbing.

Ref.	Sizes
CG1601S	S
CG1601L	L



Toes



**Little toe protective loop - CG1701**

Relieves discomfort caused by corns and calluses on the little toe.

x2		
Ref.	Sizes	
CG1701	One size	



Forefoot



**Plantar protector with pad - CG1501**

Helps reduce calluses while relaxing the metatarsal area.

x2	
Ref.	Sizes
CG1501S	S
CG1501L	L



**Extra-thin toe separator (crescent-shaped) - CG2102**

Reduce friction between the toes and abrasions that cause irritation and calluses.

x3	
Ref.	Sizes
CG2102	Mix 3 sizes



**Metatarsal pad with 2nd toe support loop - CG0701**

Helps to alleviate calluses on the sole of the foot. Reduce pressure and rubbing that causes pain when walking.

x2	
Ref.	Sizes
CG0701S	S
CG0701L	L



Forefoot



### Extra-thin non-slip mini insoles - CG0401

Provide comfort and cushioning throughout the day.

Ref.	Sizes
CG0401	One size



Sole



### Extra-thin lined insoles (metatarsal and heel support) - CCF323

Made with two densities, one especially soft for the unloading area. Prevent pressure and relieve pain on the heel and metatarsals.

Ref.	Sizes
CCF323S	S (35-37)
CCF323M	M (38-40)
CCF323L	L (41-43)
CCF323XL	XL (44-46)



### Extra-thin lined insoles (special 'all care') - CGF0370

Insoles made from gel and featuring thick shock-absorbing cushioned areas to prevent pressure on the heel and metatarsal area, relieving pain in delicate feet.

Ref.	Sizes
CGF0370S	S (35-38)
CGF0370M	M (39-42)
CGF0370L	L (43-46)

Comfort for your feet

INSOLES  
**ComforGel**

2.3

Three models to better adapt to your feet and your routines

- Everyday use
- Work
- Sport



**Every Day insoles**  
H-CGF0360 · M-CGF0361

Trimmmable. Anatomical shape and gel cushioning grooves with special width. Greater thickness in the heel area. Extra breathability. Greater comfort gel surface.

Ref.	Size
H-CGF0360	40-46
M-CGF0361	35-41



**Work insoles**  
H-CGF0364 · M-CGF0365

Trimmmable. Anatomical shape and gel cushioning grooves with special width. Flexible sides for better adaptation in stiff footwear. Breathable lining. Greater comfort gel surface.

Ref.	Sizes
H-CGF0364	40-46
M-CGF0365	35-41



**Sport insoles**  
H-CGF0362 · M-CGF0363

Trimmmable. Ergonomic cushioning with special adaptation to the foot tread area. Better foot tread stability. Breathable lining. Extra-comfort with gel grooves with special width.

Ref.	Sizes
H-CGF0362	40-46
M-CGF0363	35-41

*Gel that is not seen.  
Comfort that is felt*

ComforGel Lady insoles are presented in 4 models to better fit in all kind of women's footwear:

- Very high heels
- Daily heels
- Open shoes (or sandals)
- Flat shoes



Counter stand -  
EXPPLANTCGL

Measurements  
35 x 41,5 x 16 cm

## Comfort insoles

# ComforGel Lady



**High heels insoles -  
CG0366**

Design adapted to provide maximum comfort in the metatarsal area. Help to withstand static loads for long periods. Incorporate bubbles that provide a light massage effect when walking.

Ref.	Size
CG0366	35-40



**Daily heels insoles -  
CG0367**

Design with greater support area in the metatarsals. Prevent and relieve pain in weak, tired and delicate feet. Incorporate bubbles that provide a light massage effect when walking.

Ref.	Size
CG0367	35-40



**Open shoes insoles -  
CG0368**

Design adaptable for use with open shoes and sandals. Greater coverage in the area of the metatarsals and heel, providing maximum comfort. Incorporate bubbles that provide a light massage effect when walking.

Ref.	Size
CG0368	35-40



**Flat shoes insoles -  
CG0369**

Special design that adapts to the shape of flat shoes. Cushion and absorb impacts in the heel, knee and ankles. Incorporate bubbles that provide a light massage effect when walking.

Ref.	Size
CG0369	35-40





**Half-moon-shaped separator  
- P417**

Mild bunions or overlapping toes. Need for separation and alignment of the toes.  
Pathology: Bunions and overlapping toes.

x2

Ref.	Sizes
P417S	S
P417M	M
P417L	L



**Hourglass-shaped separator  
- P418**

Pronounced bunion with overlapping of the big and second toes, with flexibility in the big toe observed. Need for separation of these toes.  
Pathology: Bunions

x2

Ref.	Sizes
P418S	S
P418M	M
P418L	L



**Daytime bunion corrector  
- P454**

Pain from the shoe rubbing against the side of the toe. Need for protection.  
Pathology: Hallux valgus (bunions).

x1

Ref.	Sizes
P454	One size



**Night-time bunion corrector  
- P435 / P435P**

Intense pain at night and stiffness in the big toe. Need for posture correction and pain relief. Use a minimum of 2-3 times per week.  
P435: Supplied as a single unit (left or right).  
P435P: a pair (left and right).  
Pathology: Hallux valgus (bunions).

Ref.	Side
P435	L
P435	R
P435P	L and R



Heel



**Heel cup for lateral calcaneal spur - P438**

Long lasting product indicated to alleviate painful calcaneal spur and bursitis due to its two densities. Avoids the formation of callosities and provides an agreeable sensation of comfort. Central spur.  
Pathology: Calcaneal spur.



Ref.	Sizes
P438S	35 - 37
P438M	38 - 40
P438L	41 - 43



**Heel cup for central calcaneal spur - P440**

Long lasting product indicated to alleviate painful calcaneal spur and bursitis due to its two densities. Avoids the formation of callosities and provides an agreeable sensation of comfort. Central spur.  
Pathology: Calcaneal spur.



Ref.	Sizes
P440S	35 - 37
P440M	38 - 40
P440L	41 - 43



**Heel pads - P400-A**

Reduce the loads supported by the ankle, knee, hip and spine.  
Pathology: Bursitis y Tendinitis.



Ref.	Sizes
P400-A S	35 - 37
P400-A M	38 - 40
P400-A L	41 - 43



**Heel cups with orifice for calcaneal spur - P453**


Heel pain with the need to avoid putting weight on the spur area where it touches the footwear.  
Pathology: Calcaneal spur.



Ref.	Sizes
P453 S	35 - 37
P453 M	38 - 40
P453 L	41 - 43



Sole




**Non-slip anatomical insole - P441**

The metatarsal support and plantar arch place the foot in the correct anatomical position.  
 Pathology: Anatomical relief while walking. Lined to prevent slippage while walking.

x2

Ref.	Sizes
P441S	35 - 37
P441M	38 - 40
P441L	41 - 43



**Insoles with metatarsal pad - P433**

Helps relieve pain in the heel and metatarsal area. The metatarsal support and plantar arch provide comfort and relief while walking.  
 Pathology: Plantar fasciitis. Tired feet and legs.

x2

Ref.	Sizes
P433S	35 - 37
P433M	38 - 40
P433L	41 - 43



**Full gel insoles with lining - P456**

Cushion the impacts produced while walking. Lined to reduce the feeling of dampness.  
 Pathology: Tiredness in the plantar area.

x2

Ref.	Sizes
P456S	35 - 37
P456M	38 - 40
P456L	41 - 43



# 3

- 3.1. BruxiCalm
- 3.2. RespiraPlus

### Mouth guards for bruxism.

The new designs are thinner and easier to use, offering a better adaptability to the shape of mouth. More comfortable and irritation-free, for a notably improved daily wearing experience.

**BruxiCalm Day**, the guard protects your teeth against involuntary friction during times of stress (work, study, driving, etc.). • Adaptable, + anatomical • Thinner, + more comfortable • New design for improved use.



**BruxiCalm Night**, the guard protects your teeth against involuntary friction during times of stress and/or while sleeping. • Adaptable, + anatomical • Thinner, + more comfortable • New design for improved use.



**BruxiCalm Sport**, the guard protects your teeth against involuntary friction while practicing sports. • Adaptable, + anatomical • Thinner, + more comfortable • New design for improved use.



**BXCEXPOST-2**  
**Counter stand**

Ref.

827726

Measurements

32 x 44,5 x 15 cm.

Dental guard

bruxi  
calm

bruxi  
calm DAY



BruxiCalm Day, the guard protects your teeth - BXCDA

Ref.  
826036

bruxi  
calm NIGHT



BruxiCalm Night, the guard protects your teeth - BXCNIGHT

Ref.  
826040

bruxi  
calm SPORT



BruxiCalm Sport, the guard protects your teeth - BXCSPORT

Ref.  
826042

### 3.2



**Respira Plus** is a device designed to be placed in the nostrils in order to reduce resistance to incoming air and make breathing easier.

- Opens the nostrils, helping to eliminate snoring
- Improves performance in sport

Comfortable and discreet thanks to its transparency and snug fit to the nostrils.



Display box with 6 cases. -  
RP0001-6U

Ref.  
RP0001-6U

Individual case 2 pcs. -  
RP0001

Ref.  
RP0001



# 4

- 4.1. Humidifiers
- 4.2. Health care





**PRIM Umilight ultrasonic humidifier - TH001WG**

Optimises the moisture level in any room of the house, improving air quality. Essences can be incorporated.

Ideal for large spaces and/or home environments. Ambient light 5 LED colours.

Ref.	Capacity
TH001WG	2.4 l



**PRIM Bruma ultrasonic humidifier - TA-307**

Ideal for small spaces and children's rooms. Improves air quality. Essences can be incorporated.

3 in 1: Humidifier, aroma diffuser, night light (7 LED colours). Timer 1h, 2h and 6h.

Ref.	Capacity
TA-307	400 ml



**PRIM eucalyptus essence - I3132**

Essential oil for ultrasonic humidifiers. Creates an atmosphere of peace and relaxation.



**PRIM lavender essence - I3131**

Ref.	Quantity	Unid
I3132	15 m	6
I3131	15 m	6

## Health care

## 4.2



### PRIM fingertip pulse - PULSI NEGRO

- Display type: Two-colour OLED
- SpO2: Measurement range: 70-99%.
- Accuracy: 70-99%  $\pm 2$  digits 0%-69 % not defined
- Resolution: 1 %, Performance at low perfusion: PI=0,4%, SpO2. = 70% PR 30 bpm.
- Pulse: Measuring range: 30-240bpm.
- Accuracy:  $\pm 1$ bpm or  $\pm 1\%$ ., Resolution: 1bpm
- Batteries: Two AAA alkaline batteries. Consumption: less than 30mA Duration 40 hours

C.I.

820253



### PULSI INFANT PRIM fingertip pulse. Kids Oximeter - PULSI INFANT

Detailed : SpO2: Measurement range: 70% ~ 99%. Accuracy:  $\pm 2\%$  on the 80% ~ 99%,  $\pm 3\%$  on the SpO2 70%-79% Below 70% no requirement Resolution:  $\pm 1\%$ . Measurement range: 30bpm ~ 240bpm. Accuracy:  $\pm 1$ BPM or  $\pm 1\%$  (the larger one). Two OLED screen colors and different display modes. Low power consumption, allows continuous use for more than 6 hours. Rechargeable lithium battery. Battery indicator. Auto power off. Small, low weight and easy to transport. Measurements: 2.6 x 2.8 x 4.4 cm.

C.I.

821877



### PRIM arm sphygmomanometer - BU A50

Measurement of blood pressure.

For 4 different users. Stores the last 60 measurements for each user. Calculates the average for the last 3 measurements. Easily understandable colour-coded results following the WHO classification.

Ref.

BU A50



### IR thermometer - TERM\_IR

Non-contact measurement of body, ambient and fluid temperatures. Response time: <2 seconds. Measuring distance: 1 cm. 34 memories with date and time. Alarm sound in case of fever Error indicator message. Container content: 1 Infrared thermometer, 2 \* AAA batteries, 1 User Manual.

C.I.

821225



### PRIM rigid digital thermometer ECT-1

For taking under-arm and under-tongue temperatures. Measurement accuracy:  $\pm 0.1^{\circ}\text{C}$  between  $32^{\circ}\text{C}$  and  $42.9^{\circ}\text{C}$ . Previous measurement saved in the memory. 12 thermometer display stand.

Ref.

ECT-1



### PRIM flexible digital thermometer - MT519

For taking under-arm and under-tongue temperatures. Its flexible tip makes it ideal for use with children. Measurement accuracy:  $\pm 0.1^{\circ}\text{C}$  between  $32^{\circ}\text{C}$  and  $42.9^{\circ}\text{C}$ . Previous measurement saved in the memory. 12 thermometer display stand.

Ref.

MT519



**PRIM heating pad -  
WE1695PD**

Heat therapy to relieve muscle tension and joint pain.

Soft, breathable and flexible fabric. Washable: Removable control switch so that the pad and cover can be washed. Universal. Dimensions adaptable to any part of the body.

Ref.

WE1695PD



**Professional electrotherapy -  
PS2**

TENS and EMS professional electrostimulator 2-channel electrostimulator with 20 programmes. Features 6 free memories for entering desired parameters and creating individualised programmes.

Ref.

PS2



**PRIM paraffin bath -  
I3108**

Complete kit for medicinal and/or cosmetic treatment of the hands, elbows and feet.

Safe, clean and easy to use system. Includes 1 heater with intensity regulator, 2 mittens, 100 plastic bags and 2 kg of menthol paraffin wax.

Ref.

I3108



# 5

- 5.1. Mobility
- 5.2. Hygiene and daily life
- 5.3. Rest



**Rotating stand**  
12 pieces

Stand for 12 folding walking sticks

Ref.	Width	Height	Depth
600137	30cm	25cm	30cm



**Rotating combined stand**  
24 pieces

Stand for 12 folding walking sticks and 12 extendable or classic walking sticks.

Ref.	Width	Height	Depth
600283	30cm	134cm	30cm



**Wooden stand with wheels**  
17 pieces

Wooden stand for 17 crutches and/or walking sticks

Ref.	Width	Height	Depth
463811	60cm	71cm	16cm

## Walking sticks and crutches :

*Extendable anatomical walking sticks*  
Aluminium walking sticks with anatomical handle



Replacement parts



Ferrule 1,7 cm

Ref.

243860

**Extendable anatomical walking stick**

**Extendable aluminium walking stick**

Ergonomic handle

Straight bronze handle

Ref.	Side
269572	Right
533026	Left

Ref.

269572 Right  
533026 Left

600374

Characteristics



Total length

78 to 100 cm



Maximum user weight

110 kg



Material  
Aluminium

*Classic extendable walking sticks*  
Smooth aluminium walking sticks with patterns and wooden derby handle.



**Extendable aluminium walking stick**

**Extendable aluminium walking stick**

**Extendable aluminium walking stick**

Purple flower pattern

Beige flower pattern

Pink/blue flower pattern

Ref.

600070

Ref.

600071

Ref.

600075

## Walking sticks and crutches

### Classic extendable walking sticks

Smooth aluminium walking sticks with original patterns. Wooden derby handle  
It is recommended to adjust the height of the handle to the height of the user's trochanter.



**Extendable aluminium walking stick**

Black derby handle

Ref.

593814



**Extendable black aluminium walking stick**

Soft handle

Ref.

465667



**Extendable aluminium walking stick**

Bronze derby handle

Ref.

300689



**Extendable aluminium walking stick**

Turquoise pattern

Ref.

465626



**Extendable aluminium walking stick**

Beige pattern

Ref.

465610



**Extendable aluminium walking stick**

Burgundy red pattern

Ref.

465592

# Walking sticks and crutches

Luxury extendable walking sticks

Aluminium walking sticks with engraved glass effect. Wooden derby handle and matching ferrule.



**Extendable aluminium walking stick**

Blue/white flower pattern

Ref.  
465618



**Extendable aluminium walking stick**

Engraved black glass

Ref.  
600068



**Extendable aluminium walking stick**

Engraved red glass

Ref.  
600069



**Extendable aluminium walking stick**

Engraved blue glass

Ref.  
600072



**Extendable aluminium walking stick**

Engraved gold glass

Ref.  
600073

## Replacement parts



**Black ferrule**

Ref.  
243860



**Red ferrule**

Ref.  
600165



**Blue ferrule**

Ref.  
600163



**Gold ferrule**

Ref.  
600164



## Characteristics

		
Total length	Maximum user weight	Material
78 to 100 cm	110 kg	Aluminium



## Walking sticks and crutches

### Classic folding walking sticks

Smooth aluminium folding walking sticks with original patterns.

With wooden straight and derby handles. Replaceable ferrule. Include clip and bag for transportation.



**Aluminium  
folding  
walking stick**

Black derby handle

Ref.

533505



**Aluminium  
folding  
walking stick**

Bronze derby handle

Ref.

465675



**Aluminium  
folding  
walking stick**

Straight black handle

Ref.

465642

## Walking sticks and crutches :

*Classic folding walking sticks*  
 Straight handle. Include clip and bag for transportation.

*Folding walking stick with soft handle*  
 With straight padded texture handle.



**Aluminium folding walking stick**

Oriental green fantasy

Ref.

217837



**Aluminium folding walking stick**

Purple flower fantasy

Ref.

217845



**Aluminium folding walking stick**

Oriental gold fantasy

Ref.

218603



**Aluminium folding walking stick**

Black. Soft straight grey handle

Ref.

218604

Replacement parts



Round ferrules  
1.7cm

Ref.

243860

Characteristics



Total length  
85 to 95 cm



Maximum user weight  
110 kg



Material  
Aluminium

## Walking sticks and crutches

### Luxury folding walking sticks

Aluminium folding walking sticks with satin finish.

Wooden derby handle and replaceable matching ferrule. Includes clip and bag for transportation.




#### Aluminium folding walking stick

Satin pearl red flower

Ref.

600084



#### Aluminium folding walking stick

Satin pearl purple

Ref.

600085

Replacement parts



Lilac ferrule

Ref.

600162



Red ferrule

Ref.

600165



Green Ferrule

Ref.

600161

Characteristics

Total length	Maximum user weight	Material
85 to 95 cm	110 kg	Aluminium

## Walking sticks and crutches :

### Black beech walking sticks

Black beech walking sticks with derby and straight handles. Handcrafted manufacture with fine woods from sustainable Spanish forests. Replaceable ferrules.



**Black beech walking stick**  
Straight black handle

Ref.

524918



**Black beech walking stick**  
Soft straight black handle

Ref.

600080



Replacement parts



Ferrule 1,7cm

Ref.

580399

#### Characteristics



Total length

92 to 96,5 cm



Material  
Wood



**Beech mahogany walking stick**  
Straight shell handle

Ref.

527499



**Beech mahogany walking stick**  
Derby shell handle

Ref.

524900



**Beech mahogany walking stick**  
Derby honey handle

Ref.

546127



Replacement parts



Ferrule 1,7cm

Ref.

580399

#### Characteristics



Total length  
94 to 95 cm



Material  
Wood

## Walking sticks and crutches

### Aluminium walking sticks with wooden finish

Wooden finish folding or extendable walking stick with derby and straight handle. All the handles in this collection are made from Canadian maple wood.



5 height steps



**Extendable aluminium sequoia finish walking stick**  
Straight handle

Ref.  
820971



Characteristics



Total length  
71,5 a 94cm



Maximum user weight  
100 kg



Material  
Aluminium  
Puños de madera


**Extendable aluminium mahogany finish walking stick**  
Derby handle

Ref.  
820972



**Aluminium oak finish folding walking stick**

Derby handle  
Ref.  
820970



Characteristics



Total length  
84 a 94 cm



Maximum user weight  
100 kg



Material  
Aluminium  
Puños de madera

**Aluminium walnut finish folding walking stick**

Straight handle  
Ref.  
820969

## Walking sticks and crutches :



### Complete crutch - D95 6

Extendable crutch with forearm cuff and handle in impact-resistant material. Adjustable conical screw. Adjustable forearm cuff and height. Fits left and right. Pair



Ref.	Max.user weight
D95 6	125kg



### Walking stick tips - D72

Reinforce the end of the stick and cushion impacts, reducing injuries in the elbows and shoulders. Low-noise material. Pack of 10.

Ref.	Diameter
D72 3P	2.0 cm
D72 4P	2.3 cm



**Bathtub transfer board with handle - U348**

Makes bath transfer easier and more comfortable. Made from plastic. With non-slip grab handle. With holes for greater balance and drainage.



Ref.	Dimensions	Weight	Max weight
U348	73x35,5 cm	1,5 kg	136 kg



**Swivel bath chair - U276**

Ideal for patient hygiene and easy fitting. Aluminium frame and PVC chair. 360° rotation and 4 locking positions. Seat and back with drainage holes. With non-slip leg tips. Great stability. Removable armrests.



Ref.	Dimensions	Weight	Max weight
U276	70,5x41x42 cm	4 kg	109 kg



**Bath stool - U264**

Makes washing easier for the patient inside the shower without the danger of accidents. Aluminium frame with each leg height adjustable. Non-slip seat. Suction cups on leg tips for greater stability.



Ref.	High	Weight	Max weight
U264	41-51 cm	1,5 kg	120 kg



**Bath chair - U346**

Makes washing easier for the patient in a seated position without the danger of accidents. Sturdy plastic seat, EVA foam cushioning, with drainage holes for improved hygiene and comfort. With non-slip leg tips. Height adjustable.



Ref.	Height	Dimensions	Weight	Max weight
U346	36-46 cm	41x46 cm	1,6 kg	130 kg

## Hygiene - WC



**Bath chair - B1039N**

Makes washing easier for the patient in a seated position without the danger of accidents. Height-adjustable aluminium frame. Sturdy plastic seat with drainage holes for improved hygiene and comfort. Non-slip leg tips.



Ref.	Dimensions	Height	Max weight
B1039N	49x29,5 cm	36 -53 cm	100 kg



**Bath chair with back - 300686**

Makes washing easier for the patient inside the shower without the danger of accidents. Height-adjustable aluminium frame with 5 settings. Sturdy plastic seat with drainage holes for greater hygiene and comfort. Foam-lined armrests for improved grip and comfort. Easily removable back and armrests. Non-slip leg tips.



Ref.	Dimensions	Height	Max weight
300686	40x33 cm	42-55 cm	120 kg



**Bath chair with back - B1040N**

Makes washing easier for the patient inside the shower without the danger of accidents. Height-adjustable aluminium frame. Sturdy plastic seat with drainage holes for improved hygiene and comfort. Non-slip leg tips.



Ref.	Dimensions	Height	Max weight
B1040N	49x44,5 cm	36-46 cm	90 kg



**Bath chair with back - U347**

Makes washing easier for the patient inside the shower without the danger of accidents. Sturdy plastic seat, EVA foam cushioning, with drainage holes for improved hygiene and comfort. With non-slip leg tips. Height adjustable.



Ref.	Dimensions	Height	Max weight
U347	42x29 cm	35-50 cm	130 kg





**Long safety bar - U271**

Provides safety and independence. Easy to fit and very safe. Folding. Features a toilet roll holder. Lacquered in white.



Ref.	Height	Dimensions	Tube Ø	Weight	Max W.
U271	28 cm	75 cm	25 mm	2,2 kg	80 kg



**Safety bar with leg- U287**

Provides safety and independence. Folding. Height-adjustable according to needs. With 2 fastening clamps. Easy to fit and very safe. Lacquered in white.



Ref.	Height	Dimensions	Tube Ø	Weight	Max W.
U287	72-92 cm	15x75 cm	25 mm	2,8 kg	90 kg



**Clinical unisex urine bottle 'Sabot' with lid - 465850**

For people with incontinence or who are bedridden for short or long periods of time. Made from very sturdy plastic. Hygienic. Easy cleaning.

Ref.	Capacity
465850	1 L

## Others



### Large and small suction grab handles - U230 / U231

Grab handles for the safety of people with reduced mobility. Ideal for bath or shower. With safety system in case of incorrect placement. Placement on smooth non-porous surfaces.

Ref.	Size	Dimensions	Suction cup
U230	L	49 cm	8 cm
U231	S	40 cm	8 cm



### Adjustable bed backrest - U270

Ideal for bedridden people to help them read, write and eat comfortably. Hard-wearing fabric offering support and comfort. angle adjustable by means of a ratched system folding. Washable.

Ref.	Dimensions	Weight
U270	65x54x30.5 cm	2,5 kg

## Daily life



### Fridge insulin bag

Small bag/case for diabetics and their everyday paraphernalia. Includes isothermal compartment and reusable cold gel bag.

Ref.	Dimensions
265869	10x20 cm



### Stocking donner - U245

Suitable for stockings and compression tights. Made from sturdy steel tubing. Lightweight and practical.

Ref.
U245



### Plastic donner for shoes - U187

Ideal for people with difficulty bending down.

Ref.
U187



### Infant's posture pillow

Protects the head of newborn babies, up to 15 months. Designed to redistribute the pressure of the infant's head on the pillow by increasing the contact surface.

Ref.	Dimensions
246208	22x28cm



**Liquid flotation cushion - U214**

Helps prevent pressure sores on the skin and underlying tissue of people with reduced mobility. High degree of prevention and stability. Uniform distribution of pressure. Very comfortable and sturdy. Square-shaped. For chairs, armchairs, wheelchairs. Non-slip and washable outer cover. With antibacterial treatment. Crib 5 fire-retardant material.

Ref.	Dimensions	Chair width
U214 1	38 x 38 x 4,5	40 cm
U214 2	38 x 38 x 6	40 cm
U214 3	41 x 41 x 4,5	43 cm
U214 4	41 x 41 x 6	43 cm



**Anatomical viscoelastic cushion - U220**

Helps prevent pressure sores on the skin and underlying tissue of people with reduced mobility. Optimum pressure distribution. Ergonomically and anatomically shaped. With air distribution channels for better ventilation. For chairs, armchairs, wheelchairs. Non-slip and washable outer cover. With antibacterial treatment. Crib 5 fire-retardant material.

Ref.	Dimensions	Chair width	Max. user weight
U220 1	38 x 38 x 7,5	40 cm	120 kg
U220 2	41 x 41 x 8	43 cm	120 kg



**Horseshoe-shaped liquid flotation cushion - U218**

Helps prevent pressure sores on the skin and underlying tissue of people with reduced mobility. High degree of prevention and stability. Uniform distribution of pressure. Very comfortable and sturdy. Horseshoe-shaped. For chairs, armchairs, wheelchairs. Non-slip and washable outer cover. With antibacterial treatment. Crib 5 fire-retardant material.

Ref.	Dimensions	Chair width	Max. user weight
U218 3	43 x 43 x 4,5	43 cm	150 kg
U218 4	43 x 43 x 6	43 cm	150 kg



**Viscoelastic cushion density 50 with waterproof PU cover**

Relieves pressure, acting on the prevention of pressure sores in people who remain seated for long periods of time. Soft feel.

Ref.	Width	Length	Height
463084	40cm	40cm	6cm



**Wheelchair seat belt -  
574004/173899**

Prevents falls and ensures correct positioning of the patient in the chair. Pressure-fastening system for easy placement and removal. Foam velour fabric. Adjustable nylon straps.

Ref.	Dimensions
574004 M	90-165 cm
173899 L	115-200 cm



**Jacket with perineal strap -  
544262/544263**

Prevents falls and ensures correct positioning of the patient in the chair.

Ref.	Sizes	Dimensions
544262	M	100-150 cm
544263	L	120-200 cm



# 6

6.1. Plaster cast covers

6.2. Kinesiology,  
bandages and  
exercise machine

## 6.1 Plaster cast covers: Kinesiology

Safe, easy and effective method for keeping injuries dry in the shower or while taking care of personal hygiene.



### Plaster cast cover with Velcro closure

Special safety tape design

Ref.	Characteristics	Size
587014	Short arm (46cm)	Adult
300714	Long arm (70cm)	Adult
505941	Arm	Child
585182	Short leg (65cm)	Adult
584615	Long leg (90cm)	Adult
598144	Leg	Child



### Kinesiology tape

Neuromuscular kinesiology taping for the treatment of sports injuries and other physical disorders. Made from cotton and polyurethane. Hypoallergenic. Latex-free. Breathable. Waterproof.

Ref.	Measurements	Colour	Quantity
507984	5cmx5m	Red	6
507979	5cmx5m	Blue	6
507980	5cmx5m	Beige	6
507982	5cmx5m	Orange	6
507983	5cmx5m	Black	6
507986	5cmx5m	Pink	6
507985	5cmx5m	Green	6
509713	5cmx5m	Assorted colours	6

## Bandages and exercise machines

6.2



**Arm and leg exercise machine - 13029P**

Muscle rehabilitation. Suitable for the elderly or people with reduced mobility

Ref.
13029P



**Cold/heat therapy gel bag - 10210101**

Cold relieves swelling and soothes pain and heat relieves pain and increases mobility

Ref.
10210101



**Adhesive cotton bandage tape - 217870**

100% non-elastic cotton bandage tape for containment bandages.

Ref.	Measurements	Quantity
217870	3.8 cm x 10 m	3




**Beige foam pre-tape - 218579**

Soft, thin bandage tape for use under the bandage.

Ref.	Measurements	Quantity
218579	7 cm x 27 m	2







# 7

- 7.1. SPD
- 7.2. Pill organisers and crushers

# 7.1

# SPD



Pharmaceutical care tool to improve pharmacotherapeutic compliance, prevent errors and/or forgetfulness and facilitate the administration of medication.

### SPD vertical tablet

Ref.	Pieces
115840	100
115850	250

### SPD horizontal tablet

Ref.	Pieces
115824	100
115816	250

#### Characteristics



	Width	Length	Cells	Depth
115840	16cm	23cm	3x1.3 cm	2.5cm
115850	16cm	23cm	3x1.3 cm	2.5cm
115824	22.5cm	16cm	2.5x2 cm	2.5cm
115816	22.5cm	16cm	2.5x2 cm	2.5cm
115860	9.5cm	26.5cm	--	4.5cm
115865	22.5cm	16.5cm	--	3cm

### ACCESSORIES



Vertical template

Ref.  
115860



Horizontal template

Ref.  
115865

## Pill organisers and crushers



**Standard pill cutter**

Cutter that features a blade that cuts pills to obtain half the dose or enable them to be swallowed more easily.

Ref.  
599886



**Transparent deluxe pill cutter**

Completely transparent cutter.

Ref.  
492660



**Deluxe cutter and crusher**

Cuts, crushes and turns pills into powder and features 4 compartments to store them.

Ref.  
492702



**Ultra-fine crusher**

Cuts, crushes and turns pills into ultra-fine powder. Features compartments to store the pills. Special design for effortless crushing.

Ref.  
537811



**Red crusher**

Turns pills into powder without any effort to make them easier to swallow. Includes a compartment for pill storage.

Ref.  
221482

### Characteristics

	Width	Length	Depth
599886	4 cm	8.5 cm	2.5 cm
492660	3.5 cm	8.5 cm	2.5 cm

### Characteristics

	Width	Length
492702	4 cm	7.5 cm
537811	4 cm	7.5 cm
221482	5 cm	8.5 cm





# 8

- 8.1 Essential oils
- 8.2 Others

## 8.1

# OLIUM

Plants have the answer

They combine cosmetic and dermatological actions and are formulated to maintain the health and beauty of the skin, hair and other multiple uses.



SAFE DISPENSER



GLASS BOTTLE

Mixed display stand -  
OLIUMAMIX-EX12



## Cosmetic oils



### Rosehip - OLIMUM ARM-EX12

Reduces wrinkles and expression lines. Fades dark blemishes caused by exposure to the sun. Moisturises dry skin and eczema, restoring elasticity to the skin. Improves the appearance of burns and scars.

Cont.

12 PCS - 30ml



### Tea tree - OLIMUM AAT-EX9

Prevents the spread of lice and repels insects. Foot and armpit deodorant. Cleans skin with a tendency to acne and impurities.

Cont.

9 PCS - 30ml



### Argan - OLIMUM AA-EX9

Softens and restores glow to the skin. Anti-wrinkle, provides oxygenation and elasticity. Restores shine to dry and brittle hair. Provides well-being and relaxation, can be used for massages.

Cont.

9 PCS - 30ml



### Wintergreen - OLIMUM AG-EX9

Relief of muscle pain. Invigorating, combats chronic fatigue. Muscle relaxant in the bath. Fades skin blemishes, cleans pores and alleviates acne.

Cont.

9 PCS - 30ml

Measurements: 12-pc display stand: 18.8 x 22.5 x 11.7 cm

Measurements: 9-pc display stand: 14.4 x 22.5 x 11.7 cm



stickSTOP



**Alum bar - 20060101**

For closing pores and reducing razor burn during shaving.

Ref.	Pcs.
20060101	1

DETAILS



12-unit display stand (unit of sale)



**Viadol Fix Pharma (tubular mesh)**

Tubular elastic mesh for retaining all kinds of dressings, compresses or protection on any part of the body.

Ref.	Indications
70010401	Finger
70010402	Child's finger and wrist
70010403	Several adult fingers and wrist. Joints
70010404	Adult joints, hand, arm and foot. Child's leg and knee
70010405	Adult leg and knee. Child's thigh and head
70010406	Adult leg, thigh and head. Child's chest and abdomen
70010407	Adult head, chest and abdomen

## Others



**Viadol Fix Baby  
(umbilical band) - 70030201**

Tubular elastic mesh for retaining umbilical cord dressings.



**Tubular, extendable liderton**

Extendable tubular cotton bandage for protection under plaster.

Ref.		Width	Length
516534	No. 1	1.2cm	20m
504688	No. 2	2cm	20m
504779	No. 3	2.8cm	20m
490045	No. 4	5cm	20m
508549	No. 5	7.5cm	20m
597286	No. 6	10cm	20m
529693	No. 7	15cm	20m
146381	No. 8	20cm	20m



**Tubular, compressive lidergrip**

Compressive tubular elastic bandage for rehabilitation support.

Ref.		Width	Length
515585	A	4.5cm	10 m
525519	B	6.3cm	10 m
542894	C	6.8cm	10 m
573816	D	7.5cm	10 m
542282	E	8.8cm	10 m
515882	F	10cm	10 m
572545	G	12cm	10 m
146407	J	17.5cm	10 m
564849	K	21.5cm	10 m
146399	L	32.5cm	10 m



**PRIM**

**PHARMA**

general catalogue

**PRIM S.A.**

Polígono Industrial nº 1  
Calle F, 15  
28938 Móstoles. Madrid

<https://international.prim.es/>