




Finger extensor splint^{PSL180}

Dynamic splint for proximal interphalangeal joint extension.


Helps immobilise and reduce stress on the joint, alleviating the flexor tendon injury.



Flexor tendon injuries occur very frequently in accidents at home, at work or in sports.

These injuries may lead to complications and early rehabilitation and **immobilisation** are essential to avoid them.

Dynamic splints relieve pain and favour proper healing, which helps prevent further injury.





Finger extensor splint PSL180

Ergonomic comfort



adaptive stabilising ring

Stabilises the IFP joint, preventin deviation and unwanted movement.

3-point ergonomic support

enhances the user's comfort while keeping the finger in extension.

dynamic vertical structure

Enables flexibility, facilitating joint movement while rehabilitation is carried out.

Option to increase or reduce the tension depending on the degree of injury.

Padded for total dermatological compatibility.




* adaptable to the user's anatomy and inflammation phase

Size	Sizes	
	A Range (cm)	B Range (cm)
XS	5 - 6,3	1,3 - 2,3
S	6,3 - 7,6	1,3 - 2,3
M	7,6 - 8,9	1,3 - 2,3
L	8,9 - 10,2	1,3 - 2,3



Range A: from the metacarpal head to the distal interphalangeal joint



Bilateral

100% LATEX FREE ✓