


Herbitas
Laboratorios

STEP TO THE FUTURE


newfeet

**PRINTING SYSTEM
OF 3D INSOLES DualPRO**






DualPRO



- *Newfeet DualPRO improves printing capability compared to previous systems.*
- *It prints vertically, enabling molded base insoles.*



REPLICATE A INSOLE WITH COMPLETE ACCURACY



**INCREASED PRODUCTION VOLUME I
N LESS TIME**




AUTONOMOUS MANUFACTURING



DUST-FREE WORK AREA



ZERO WASTE



2 year warranty




- ✓ Simpler and more intuitive software
 - ✓ New file with customer data
- ✓ Vectorized template for greater precision in the final design

SOFTWARE VERSIONS AVAILABLE

○ **Standard**

○ **Premium**




Improved software that allows the creation of templates automatically with the maximum guarantee of adaptation.

NEW FEATURE

1 *Take a photo of the insole base*

2 *The software generates the last from the captured image*

3 *You can use this last as many times as you like*





Standard scanner



Portable scanner



*Real image scanning

Scanner with textured image quality

The Newfeet scanner incorporates a processor that allows generating images with greater speed and quality.

Scanned with realistic images like you have never seen, which allows for better diagnosis.




DualPRO



WIFI CONNECTION

NEWFEET SPECIAL FILAMENT

We are manufacturers. Our filament offers special features that guarantee exceptional durability, flexibility, and comfort, perfectly adapting to the needs of orthopaedic and sports footwear. With our advanced technology and high-quality materials, we ensure that each print is precise and meets the most rigorous demands of design and functionality



Multiple densities



High durability



Allows for lining



Flexibility



Excellent memory



Breathability



Polishable



Lightweight

AVAILABLE COLORS FOR NEWFEET FILAMENTS

STANDART



BLACK



GREY



ORANGE



CRYSTAL



GREEN



YELLOW



WHITE



PINK



LILA

SOFT



BLUE SOFT



CAMEL SOFT

RIGID



CARBON



POLYPROPYLENE

The logo for newfeet features a stylized white graphic of three wavy lines above the brand name "newfeet" in a lowercase, sans-serif font.

DualPRO



Herbitas
Laboratorios

STEP TO THE FUTURE

Laboratorios Herbitas S.L.
Periodista Badía, 13B
46134 Foios (Valencia)
España

Tlf: +034 963627900

herbitas@herbitas.com

herbitas.com