



NEW Post-surgical shoes

Comfortable and safe protection

Post-surgical shoes that protect and keep the foot in the correct position during the recovery process.





Comfortable and safe protection

NEW

Post-surgical shoes

PRIM's new post-surgical shoes feature a Velcro closure system that can adjust to feet with significant deformities or post-operative bandages or plaster casts.

Made from breathable fabric to ensure comfortable use.



non-slip
sole



fits left
and right



extensive
size range

• PS100

NEW



• PS200 open with off-loading wedge

With extra ankle support
Ergonomic sole

NEW



• 9115 closed with off-loading wedge



PS100 PS200



shoe number

sizes	(0) XS	(1) S	(2) M	(3) L	(4) XL	(5) XXL
	33/35	36/38	39/40	41/42	43/45	46/48
	● ●	● ●	● ●	● ●	● ●	●

9115

sizes	(0) XS	(1) S	(2) M	(3) L	(4) XL	(5) XXL
	35/36	37/38	39/40	41/42	43/44	45/46
	●	●	●	●	●	●

