



Ecoprim *continuous*
compression **thoracic protector**

100% LATEX FREE ✓




Thoracic
protector

Ecoprim

Both models are indicated for patients recovering from cardiac or thoracic surgery.

This product is also recommended for patients with neurological diseases, neurodegenerative diseases, dementia, etc., where the patient is unable to manage compression on demand independently.


They ensure adequate protection for patients who are unable to apply compression on demand during their recovery process..



Made in
open-cell
mesh fabric



Breathable



Adaptable to
the patient's
anatomy



Elasticated
and
non-elastic
materials for
optimal fit

Ecoprim continuous compression thoracic protector *female*

ECOPRIMCCM

Ecoprim continuous compression **thoracic protector** specially designed for women.

The central **pad** protects the wound

The **straps** pass over and under the breasts, **ensuring a precise fit** and allowing **correct regulation** of the compression level

Push-and-click buckles for rapid opening

- Provides **compression to the thorax** without affecting the breasts so that patients can lead a normal life during the recovery period.
- The **central pad** protects the wound from any unwanted pressure or chafing.
- The straps go around the breasts, leaving them free from **pressure and constraint at all times while ensuring the necessary compression is applied to the chest.**
- **Intuitive waistcoat design** – easy for the patient to put on and take off.
- **Padded sections** simultaneously cushion and compress for **enhanced comfort.**
- The **central buckles** are positioned over the central pad to avoid any possible contact with the wound.



Ecoprim *continuous compression* *thoracic protector unisex*

ECOPRIMCCU

Ecoprim unisex
continuous compression
thoracic protector.

Side pads that rest on the
pectoral muscles

Straps that ensure a precise fit
and allow correct regulation of
the compression level

Push-and-click buckle for
rapid opening

Provides **compression to the thorax**, allowing patients to lead a normal life during the recovery period.

Can be used by both men and women.

The **two side pads** sit on the pectoral muscles, leaving the area around the wound free.

Intuitive waistcoat design – easy for the patient to put on and take off.

Padded sections simultaneously **cushion and compress** for enhanced comfort.


Fastens with **two central buckles** positioned over the side pads, preventing them from coming into contact with the patient or exerting any form of pressure on the wound.



These **new protectors** are designed **so that the specialist can adjust the compression level before the patient goes home**. Once at home, the patient can put it on and take it off without changing the compression level set by the specialist.

The **new design** eliminates the need for sudden or **awkward movements** when putting on or taking off the protector. The shoulder straps fasten with buckles for an effective, secure and comfortable fit.

The thoracic protector has **push-and-click buckles positioned over the lateral chest area** to prevent **pressure from being applied to the wound**. In an emergency, the protector can be opened in seconds by squeezing the sides of the buckles.



Sizes and measurements

ECOPRIMCCM *female*

Size	Chest contour measurement
1	70 - 100 cm
2	100 - 130 cm


Fits left or right

Both models are available in two sizes

ECOPRIMCCU *uniSex*

Size	Chest contour measurement
1	75 - 105 cm
2	105 - 135 cm


Fits left or right





C/ Yolanda González, 15 - Pol. Ind. N°. 1
28938 Móstoles - Madrid - SPAIN

Tel. : (+34) 913 342 532 / (+34) 913 342 539
e-mail: export@prim.es



<https://international.prim.es/>